ELEKTROCHIRURGIE ELECTROSURGERY

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Electrosurgery is a technique used for nearly all types of surgical procedures nowadays. In order to use electrosurgery effectively, it is paramount to understand it and to learn and apply the principles that guarantee maximum safety for both the surgeon and the patient.

SURGICAL DIATHERMY

A surgical diathermy generator is a device that transforms electrical energy into high frequency alternating current (HF).

Flowing through the tissue, HF current, produces a thermal effect used for cutting and coagulating tissues.

A surgical diathermy generator produces alternating current of the frequency above 300 kHz, which prevents the risk of accidental electrolysis or muscle and nerve stimulation.

MONOPOLAR OPERATION

In the monopolar mode, HF current is applied to the tissue with an active electrode. The effect of section or coagulation is possible thanks to the concentration of a high density of HF current on a small surface of the active electrode. This results in a temperature rise and the evaporating of water from the tissue directly surrounding the active electrode, thus creating the effect of haemostasis and cessation of bleeding or cutting.

HF current then flows to the neutral electrode and is dispersed on its surface. This enables lowering the density of HF current, and there is no undesired thermal effect where the neutral electrode is applied. Through the neutral electrode, the current flows back to the generator.

BIPOLAR OPERATION

During bipolar operation, HF current flows between the two sides of a bipolar instrument and is only concentrated on the small surface of tissue between them.

Bipolar work carries no danger of the current flowing through the patient's body. The risk of burns in areas distant to the operation is minimized.

Bipolar coagulation is therefore safer than monopolar operation and is particularly recommended for patients with a pacemaker or for operations on organs with a small cross-section. There is no need to use the neutral electrode for bipolar operation.


While using surgical diathermy generators producing HF current, it is necessary to remember two basic principles:

- The current flows through all accessible routes.
- HF leakage current flows between two nearby conductors, even if these do not touch each other.

■ PREPARATION FOR THE PROCEDURE AND POSITIONING THE PATIENT

- The patient should lie on a dry and insulated surface. The patient may not touch any metal parts. If during the procedure the surface gets wet, it should be exchanged for a dry one, as long as it is feasible. No liquids should be accumulated beneath the patient.
- Parts of the patient's body that touch each other, such as fingers, limbs, should be insulated. This can be done by separating them with dry gauze or surgical paper. Such insulation will prevent potential burns.

- Surgical diathermy generators may interfere with the operation of other electric devices present in the theatre, for example EEG monitors. These should be placed at a distance from the generator. HF cables should not be led together with cables powering other electronic devices. The electrodes of the monitoring systems should be placed at least 20 cm from the operation area and the area of neutral electrode application.

- Only non-flammable disinfectants may be used. Should the need arise to use alcohol-based solutions, time should be allowed for them to evaporate, before the procedure. There is also the risk of spilling such solutions beneath the patient's body or onto it. Should this occur, the spillages must be thoroughly dried. During electrosurgical procedures, the sparks created at the active electrode create the risk of igniting the dressing materials, as well as metabolic gases.

APPLICATION OF THE NEUTRAL ELECTRODE

- Before application, the state of the neutral electrode and the connecting cable should be inspected.
Do not use if there are any visible signs of defect or damage to the insulation.

- The neutral electrode should be applied in a way that allows it to adhere completely to the patient's body.

- A reusable silicone electrode should be fastened in a way that makes it impossible for it to relocate - a special tape should be used to fasten the neutral electrode.
- For procedures involving small children, special sizes of electrodes should be used: for babies and children.

- A neutral electrode must not be splashed with any liquids or wrapped with anything.
- No extra electrically conducting gel is to be put on top of the neutral electrode.
- Remember that no liquids can remain between the electrode and the patient's body.

- While disconnecting the neutral electrode, the cable must never be pulled.
- Under no circumstance may an attempt be made at repairing a neutral electrode.
- A neutral electrode should be placed close to the area of operation, but at least 20 cm away from it.
- The neutral electrode should be placed in such a way that the longer edge is directed at the area of operation - this does not apply to the disposable TwinSafe electrodes, which can be applied in any direction.

- Before applying the neutral electrode, familiarize yourself with the instruction manual provided by the manufacturer.
- The manufacturer's instructions on the package of the neutral electrode must always be adhered to.

- All **HEBU** generators are equipped with a monitoring system for the neutral electrode **NEM** (Neutral Electrode Monitor). Using **disposable split neutral electrodes** combined with the NEM system guarantees constant control of the electrode's adhesion and ensures maximum security of the procedure.

- With disposable electrodes, the expiry date must always be checked. Electrodes that are past expiry date must not be used.
- A disposable electrode may only be used once.
- Disposable electrodes must be applied with care. Once applied, an electrode should not be reapplied onto a different place. If it is necessary to change the place of application, a new electrode should be used.
- While lifting-off disposable electrodes, attention must be paid not to damage the patient's skin. No sudden movements should be made, the electrode should be lifted gently.


- After each change of the patient's position, the application of the electrode, and the connecting cables must be checked.

PLACES OF APPLICATION OF A NEUTRAL ELECTRODE

- A neutral electrode must always be applied onto clean and disinfected skin.
- Neutral electrodes should be applied on places with high vasculature, with no hollows or bending of the skin occurs, for example an arm or a thigh.

- Do not apply the electrode onto scars, cuts or scratches.
- Do not apply on protruding, bony surfaces.
- Do not apply the electrode on surfaces with a lot of hair-if necessary, shave the area of application.
- Do not apply on extensively fatty areas, such as the stomach or the rear.
- Do not apply on surface of skin under which there are implants.

CABLES AND ACTIVE ELECTRODES

- Before beginning the procedure, the state of the cable insulation and active tips should be inspected. **DO NOT USE DEFECTIVE CABLES AND INSTRUMENTS.**
- Do not attempt to repair instruments or connecting cables. Substitute damaged electrodes with new ones.
- Before beginning the procedure, check the connections of the cables and the application electrodes. It is necessary to ensure all the equipment works properly.


- The connecting cables should not be rolled up, fastened or attached to the operating table with metal clamps. The cables should be led in a way that avoids touching the patient and other cables.
- The handle with the active electrode must never be put directly on the patient's body. After performing the procedure, put it aside, this will prevent accidental activation of the generator, which can result in burns.
- The generator should not be activated if the active tip is not placed on the spot of the intended application. This can result in damaging other tissues.

- Special care must be taken while using the active electrode in areas where flammable substances may be present, e.g. disinfectants or alcohol based anaesthetics.
- For patients who are connected to monitoring devices (EKG), the monitoring electrodes should be placed at the maximum possible distance from the areas of application of the electrosurgical electrodes.


- Dirty active electrodes may adversely influence the quality of work of the generator.

- It is recommended to clean the active electrodes during the procedure - all the tissue residue on the surface of the electrode decreases the effectiveness of the operation. Disposable cleansers are recommended for cleaning electrodes.
- To disconnect the cables and neutral electrodes, always grab the plug. **DO NOT PULL THE CABLES.**

GUIDELINES ON USING ELECTROSURGERY

- Before using a surgical diathermy generator, we recommend familiarising the staff with the effects of different operation modes and settings of power, by performing a trial on fresh, fibrous meat (beef, pork).
- The lowest effective power setting should always be used for a given procedure.
- During the procedure, there is a small risk of interfering with the operation or damaging the pacemaker. In such cases, bipolar technique is recommended. Before using electrosurgery, consult an authorised representative of the pacemaker manufacturer.
- **Using electrosurgical generators for out-patients with pacemakers is not acceptable.**
- While using spray coagulation, maintain safe distance from fingers, metal parts of the endoscope and fibrescope optics.


- The often applied technique of coagulation “via tools” can only be applied with the use of appropriately insulated pincers. There are special pincers with insulated handles (more information can be obtained from **HEBU** representatives). Surgical gloves do not provide the operator with sufficient protection from burns. **Never use spray coagulation with this technique.**

- In order to prevent the tissue from adhering to the active electrode (mono or bipolar), it is recommended to wait a moment after coagulating the tissue and turning off the current flow, to allow the active tip to cool down, before removing the electrode from the tissue.
- Remember that while coagulating vessels or organs with a relatively small cross-section, the density of the current increases and may lead to burns - in such cases bipolar coagulation is recommended.
- **Do not use bipolar scissors for the bipolar cutting mode. When using the scissors, apply bipolar coagulation.**

ARGON COAGULATION

- When using argon coagulation, take all the precautions applicable for standard monopolar coagulation.
- Familiarise yourself with the instruction manuals for argon accessories.

- Before beginning the procedure, the condition of the argon applicators should be thoroughly inspected. Do not use applicator if the insulation is broken.
- Before activating, argon instruments must be filled with gas, by using the mode **PURGE** of the argon diathermy generator.
- Do not place the argon applicator directly on the tissue. The distance between the applicator tip and the tissue should be 3 to 5 mm.
- While activating, do not touch vessel membranes of adjoining tissues with the applicator tip.
- Do not blow argon into the vessel system.

ARGON IN LAPAROSCOPY

- During laparoscopic procedures the argon flow causes the rise of the insufflation pressure. PERFORMING ARGON LAPAROSCOPIC PROCEDURES IS ONLY POSSIBLE WITH INSUFFLATORS THAT HAVE A PRESSURE LEVELLING FUNCTION. In case of any doubts, please contact the insufflator supplier in order to confirm that such a function exists.

- In order to prevent a sudden rise of insufflation pressure during the application of argon, the trocar valve should be open. If the pressure attains critical level, stop the application of argon and wait until it drops below the critical level.
- Apart from monitoring the pressure of the pneumoperitoneum with the use of the insufflator, an independent, constant control of the pressure is required on the part of the operating team.

- Argon 4.5 (99.999%) is used for argon coagulation

CLEANING AND DISINFECTION

- Remember that the reusable instruments attached to the generator are not sterile.
- All the electrosurgical accessories should be washed, sterilized and stored strictly following the manufacturer's recommendations.
- Remember that all the reusable accessories have a limited life and they are subject to wear as they are applied.

STERILIZING THE ACCESSORIES

- Sterilization methods should follow the recommendations of the supplier of given accessories. The accessories supplied by the manufacturer, unless otherwise marked, are not sterile and require sterilization before use.
- The accessories offered by **HEBU**, unless marked as disposable, may be sterilized in steam, with the maximum temperature of 134°C and the pressure of 2.3 Bar or in gas.

SILICONE NEUTRAL ELECTRODES

- Return electrodes should be washed and disinfected before use.
The disinfectants available for return electrodes: Incidin perfekt, Incidur F.
- With reusable neutral electrodes, it must be remembered that due to the washing out of the active substances while cleaning the electrode, silicone loses its conductive properties with time. In order to ensure the security of procedures, reusable electrodes should be serviced periodically together with the generator.

CLEANING THE GENERATOR

- The generator should be kept clean. It is easy to keep the device clean, it may be cleaned with standard disinfectants. Avoid surplus cleaning liquids penetrating inside the device.

SECURITY SYSTEMS FOR HEBU ES GENERATORS

- All HEBU electrosurgical generators are marked with the CE label.


■ AUTOTEST

After turning the generator on, an internal test is performed that includes all the subassemblies of the generator and the accessories connected to it.

■ NEUTRAL ELECTRODE MONITORING SYSTEM - NEM

HEBU generators are equipped with an application monitoring system for disposable split electrodes. Constant monitoring of the application during the procedure guarantees maximum security. If the application of the electrode is unsatisfactory, the generator stops further operation.

For reusable silicone electrodes, appropriate connection of the electrode to the electrosurgical generator is monitored.

■ AUTOMATIC POWER ADJUSTMENT

HEBU generators are equipped with an automatic power adjustment system. Thanks to microprocessors that constantly supervise the processes of cutting and coagulation, the generator automatically adjusts the power within a given range, in order to achieve the optimal effect.

■ OVERLOAD

HEBU generators are secured against excessive activation time.


Generator overload is signalled by a red OVERLOAD control light. The generator then induces a cessation of the procedure until the control light goes off.

■ SECURITY AGAINST DEFIBRILLATION IMPULSE

HEBU generators are manufactured in class I CF and are secured against defibrillation impulse, according to the **EN 60601-1**.

■ SERVICE CODES

In case of an abnormality, the generator displays a message - a service code. The index of the codes can be found in the instruction manual supplied together with the generator.

- The manufacturing process of HEBU generators is based on state-of-the-art technologies and maintained according to the highest standards, which is confirmed by the certifications of the Quality Management System **EN ISO 9001:2000**, **EN ISO 13485:2003** and **93/42/EWG**.

SERVICE AND PERIODIC REVIEWS

- With the view of increasing security of electrosurgical procedures, the installation of an **HEBU** generator is always combined with a **free user training session** covering the areas of using the electrosurgical generator and the principles of security in electrosurgery.
- The training for users is performed by an **HEBU** representative and ends with a certificate.

■ **HEBU generators have a 24-month guarantee**

- The users of medical equipment are obliged to take due care while selecting, installing, activating, performing service and maintenance, and in particular to follow the instruction manuals supplied by the manufacturer.
- **HEBU** as a manufacturer of electrosurgical generators recommends conducting periodic reviews **at least once a year.**

- Caring about the safety of our Customers, we inform that periodic reviews, as well as guarantee and post-guarantee repairs of **HEBU** electrosurgical generators may only be performed by units authorised by **HEBU**. Service authorisation confirms professional technical skills needed to perform such a service.

- We provide professional technical and substantial help to our Customers.

HEBU medical GmbH

Badstraße 8
78532 Tuttlingen / Germany

Tel. +49 74 61 94 71 0

Fax +49 74 61 94 71 22

info@HEBUmedical.de

www.HEBUmedical.de