HF 9530-02

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HF 9530-02

Contents

1	Used symbols.....	11
2	Warning.....	12
3	Introduction	12
4	Scope.....	12
4.1	Inteded Use	12
4.2	Contraindication.....	12
5	Product description.....	13
6	Commissioning and use	14
6.1	Connection to gas cylinder.....	14
6.2	Connection to HF-Device.....	14
6.3	Use.....	15
6.3.1	Opening.....	15
6.3.2	Closing.....	15
6.3.3	Adjusting the mechanical style pressure switch adjustment screw	15
7	Cleaning.....	15
8	Maintenance.....	16
8.1	Changing gasket ring.....	16
8.2	Changing filter screw	16
9	Storage.....	16
10	Technical data.....	16
10.1	Pressure reducer	16
10.2	Mechanical style pressure switch.....	16
11	Warranty / Repair	17
12	Manufacturer and servicing address.....	17

1 Used symbols

Please read these instructions carefully and keep them for future reference. Improper handling can cause personal injury or damage the product.

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Additional symbols used:

Symbol	Definition
	CE marking
	Attention
	Manufacturer

2 Warning

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

!	Incorrect use of the pressure reducer can cause serious damage.
!	Avoid mechanical overstressing of the device beyond the structural design, this can lead to breakage and loss of function!
!	Run the pressure reducer only at the intended pressure and with the intended gas. The use of improper gas or at improper pressure might cause risks of explosion.
!	Protect the pressure reducer from dust, oil and other sources of dirt.
!	Do not use the pressure reducer if it is not in perfect working conditions.
!	When you draw gas, the cylinder must be placed upright and protected from falling.
!	Danger of explosion. At startup and also during use, hands and tools must be free from oil and grease.

3 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

4 Scope

This product is designed to power HBS-Xtouch devices with argon gas. It reduces the gas pressure of the gas source to the required operating pressure and maintains it constant.

4.1 Intended Use

Argon pressure reducer: This product is designed to power Argon gas HBS x-touch devices. It reduces the incoming gas pressure of the gas source to the required operating pressure and maintains it constant.

4.2 Contraindication

No contraindications are known.

5 Product description

- High pressure gauge (1)
- Mechanical style pressure switch (2)
- Pressure adjusting screw (3)
- Relief valve (4)
- Outlet connection (5)
- Body (6)
- Gasket ring (7)
- Coupling nut (8)
- Inlet connection (9)
- Filter screw (10)
- Adjusting screw mechanical style pressure switch (11)


Figure 1

Figure 2

6 Commissioning and use

6.1 Connection to gas cylinder

The pressure reducer is constructed for gas cylinder connections in line with DIN 477.

- At every gas cylinder change, initial check the gasket ring (7) on intactness.
- Before connecting the pressure reducer to the gas cylinder, make sure the openings of the gas cylinder and the pressure reducer are free from dirt. Otherwise blow out with compressed-air.
- Turn the pressure adjusting screw (6) anti-clockwise as far as it will go. Do the same with the valve of the gas cylinder.
- Put the inlet connection (9) of the pressure reducer on the opening of the gas cylinder. Align the high pressure gauge (1) upright and close the coupling nut (8) manually. If necessary again align the high pressure gauge.
- Tighten the coupling nut by using a 30 mm flat wrench.

Stay away from heat sources and easily inflammable materials during the commissioning.

Do not direct the opening of the gas cylinder on yourself or other persons.

6.2 Connection to HF-Device

Ports of the HF-Device


Figure 3

- gas pressure measure line 2 (12)
- gas connection 2 (13)
- gas pressure measure line 1 (14)
- gas connection 1 (15)

- Plug in the mechanical style pressure switch (3) leaving cable in gas measure line 1 or 2 (14 or 12).
- Turn the lock clockwise until it is engaged. You can hear that at the click-sound.
- Plug in the outlet connection (5) leaving tube with its adapter at the end to gas connection 1 or 2 (15 or 13). The lock engages by itself (Figure 4 and 5).


Figure 4


Figure 5

6.3 Use

6.3.1 Opening

- Slowly open the gas cylinder valve. The high pressure gauge (1) will show you the cylinder pressure.
- Slowly open the pressure adjusting screw (3) as far as it will go.

Opening the gas cylinder valve too quickly might cause gauges malfunction.

6.3.2 Closing

- Close the gas cylinder valve.
- Release the gas out of the pressure reducer.
- Turn the pressure adjusting screw (3) anti-clockwise till it is completely closed.

6.3.3 Adjusting the mechanical style pressure switch adjustment screw

The mechanical style pressure switch fulfills the purpose to provide the control system of the HF-device with information, whether the required pressure is on or not. For that you can set a minimum pressure with the aid of the adjustment screw. Pay attention that the minimum pressure is always set under 4 bar, otherwise the gas supply will be inhibited.

To set the minimum pressure, perform the following instructions:

- Turn the adjustment screw (11) clockwise with the aid of a slotted screwdriver, to increase the pressure.
- Turn the adjustment screw (11) anti clockwise with the aid of a slotted screwdriver, to decrease the pressure.
- The adjusted pressure can be seen on the scale.

7 Cleaning

In case you need to clean the pressure reducer use a dry, dust-free cloth to wipe the item. For stubborn soiling use a damp cloth.

Avoid water inside the pressure reducer. It might cause risks of explosion.

8 Maintenance

Hands and tools must be free from oil and grease during maintenance

8.1 Changing gasket ring

At a high level of abrasion perform the following instructions:

- Lever the gasket ring with the aid of a slotted screwdriver and take it off.
- Pull the fresh gasket ring without any aid of tools on its place.

8.2 Changing filter screw

At a high level of soiling, but at least every 2nd year perform the following instructions:

- Release the filter screw with the aid of a 6 mm allen wrench.
- Hold the pressure reducer so that the inlet connection (9) is facing to the ground. This way it is prevented that dirt enters the inside.
- Turn out manually the filter screw.
- Turn in manually the fresh filter screw.
- Tighten the filter screw with the aid of a 6 mm allen wrench.

9 Storage

The pressure reducer must be stored dry and dust-free. The surrounding air must be non condensing.

It is recommended to keep the wrapping, so the item can be storage in it when required.

10 Technical data

10.1 Pressure reducer

Inlet pressure P1 [bar]	max. 230
Outlet pressure P2 [bar]	4
Delivery flow Q1 [m ³ /h]	< 1
Operating temperature [°C]	-20 bis +60

10.2 Mechanical style pressure switch

Model	IS 1000		
Proof pressure [bar]	10		
Max. pressure [bar]	7		
Operating pressure range [bar]	1 bis 6		
Hysteresis [bar]	0,8		
Ambient and gas temperature [°C]	-5 bis +60		
Max. contact capacity	2 VAC AC / 2 W DC		
Voltage	24 V AC/DC	48 V AC/DC	100 V AC/DC
Max. operating current	50 mA	40 mA	20 mA

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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