*HF 9580-00*  
*HF 9580-34*  
*HF 9580-01*  
*HF 9580-02*  
*HF 9580-03*  
*HF 9580-04*  
*HF 9580-05*  
*HF 9580-99*  
*HF 9928-13*

**Bipolare Koagulationszange für die Laparoskopie /**  
**Bipolar coagulation forceps for laparoscopy /**  
**Bipolar coagulation forceps for laparoscopy**

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Seite 3

**Operating Manual**

Page 9

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Page 15

**DEUTSCH**

**ENGLISH**

**ITALIANO**

*HF 9580-00*  
*HF 9580-34*  
*HF 9580-01*  
*HF 9580-02*  
*HF 9580-03*  
*HF 9580-04*  
*HF 9580-05*  
*HF 9580-99*

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# Bipolar coagulation forceps for laparoscopy- HEBU HF 9580

## Contents

1	Application and Safety.....	10
2	Scope .....	11
2.1	Intended Use .....	11
2.2	Contraindications.....	11
3	Assembling and Disassembling.....	11
3.1	Disassembling .....	11
3.2	Assambling .....	12
4	Reprocessing .....	13
4.1	Mechanical Reprocessing .....	13
4.2	Manual Reprocessing.....	13
5	Storage.....	14
6	Warranty.....	14
7	Manufacturer and Service address .....	14

# 1 Application and Safety

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

!	All instruments must be before first putting into use before any further Application completely cleaned, disinfected and sterilized.
!	It is very important to each surgical instrument before each use visible Damage and wear, for example Cracks, fractures or failures to examine the insulation. In particular areas such as blades, points, locks, ratchets and latches, and any moving parts, insulation and ceramic items are carefully scrutinized.
!	Do not use damaged instruments
!	Do not use in the presence of flammable or explosive substances.
!	The instrument should not be placed on the patient.
!	Coagulate only when the contact areas are within view. Don't touch any other metal instruments.
!	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
!	Avoid the simultaneous use of HF surgery and electrically grounded electrodes and sensors of diagnostic devices!
!	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
!	Observe manufacturer's instructions for electrical load capacity!

## 2 Scope

The use of HF bipolar clamps is used in HF surgery for coagulation and tissue cutting. In addition, for the permanent closure of veins, arteries and tissue bundles with subsequent separation.

### 2.1 Intended Use

HF coagulation forceps, laparoscopy; HF Clamps: A bipolar product covered with electrical insulation down to the working end. It is used in open surgery to coagulate vascular or tissue bundles for hemostasis. It is connected either to an electrosurgical diathermy unit or to an electrocautery generator that delivers the bipolar electrical current. It is typically designed as a scissors-type instrument with a long, curved, pincer-like mouth. It is made of metal and synthetic material in different sizes and shapes. It is a reusable product.

The Instrument is inserted through a 5.5 mm trocar sleeve. Coagulation results out of electrical energy generated by RF generator for electrosurgery. This instrument can be connected to the bipolar output of RF generators of the HEBU HBS x-touch series. Using a suitable bipolar it is connected to the output of a bipolar RF generator. If indexed, specifically bipolar Coagulation can be used. The maximum output voltage of the generator must not exceed 500 V<sub>pp</sub>.

### 2.2 Contraindications

The RF bipolar clamps should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.

## 3 Assembling and Disassembling

### 3.1 Disassembling


Figure 1


Figure 2

### 3.2 Assambling


Figure 3


Figure 4

## 4 Reprocessing

### 4.1 Mechanical Reprocessing

Because of product design, materials used and the intended use is no defined limit of a maximum feasible Treatment cycles should be established. The life of the instruments is determined by their function and the careful use of them determined. Instruments for Electro-Surgery under are naturally increased wear, depending on the type and duration of use.

#### Preparation

Immediately after application remove coarse dirt from the instruments. No fixation agents or hot water (> 40 ° C) using, since the fixation of residues and cleaning can affect success.

**!!! Also important thread of the inside of the mouth part!!!**

The handles should be thoroughly rinsed with the flush adapter (HF 9580-99).

Do not use any kind of acidly cleaning agents.

#### Disinfection:

Automatic disinfection, taking into account the national requirements for the A0 value (see ISO 15883).

#### Drying:

Dry the outside of the instruments by the Drying cycle of the washer disinfectors. If necessary, in addition, a manual drying with a lint-free cloth can be achieved. Cavities of instruments with sterile Air dry.

### 4.2 Manual Reprocessing

#### Cleaning:

**!!! Instruments must be disassembled before cleaning!!!**

Prepare a cleaning under Manufacturer details.

- Rinse the instruments under cold tap water (<40 ° C) until all visible Contaminants are removed. Imbedded dirt is removed with a soft brush
- Insert the instruments completely in the prepared cleaning. Exposure under manufacturer's instruction.
- Clean the instrument inserted manually with a soft brush. All surfaces must be brushed several times.
- Applies only for channels and internal surfaces: Take the brush at least six times in the tubes inside and out. Rinse then the tubes with distilled water and repeat this procedure.
- Rinse the products very thoroughly under running tap water for complete removal of detergent.

#### Disinfection:

- Prepare a disinfectant according to information provided by the disinfectant.
- Place the instruments in the disinfectant, and observe the prescribed contact time.
- Rinse the products very thoroughly with deionized water to complete removal of the disinfectant.

**Drying:**

The manual drying is effected by a lint-free cloth and, in particular for drying cavities and canals with sterile compressed air.

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Visual inspection for cleanliness, if applicable, installation and function test according to the manual. If necessary, repeat the treatment process until the instrument is visually clean. The original shine can be restored by rubbing with a cloth.

**Steam sterilization:**

Sterilization of instruments by fractional pre-vacuum procedure (according to ISO 13060 / ISO 17665) take into account the national requirements.

- 3 pre-vacuum with at least 60 mbar
- Heating to a sterilization temperature of at least 132 ° C, max. 137 ° C
- Shortest time: 3 min.
- Drying time: min at least 10.

## 5 Storage

Storage sterilized instruments in a clean, dry and dust-free environment at moderate temperatures of 5°C to 40°C.

## 6 Warranty

We provide a guarantee for all HEBU products in accordance with statutory and country-specific legislation (proof of purchase by invoice or delivery note). The warranty period is two years.

Damage caused by incorrect handling or normal wear and tear is not covered by the warranty.

In case of warranty claims, please enclose a copy of the invoice or delivery note with the device.

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**HEBUmedical GmbH**