HF 9581-20
HF 9581-37

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HF 9581-20

HF 9581-37

Vessel sealing instrument BiSeal[®] with cutting function

Table of Contents

1	Directions for user manual	13
2	Application and Safety	13
3	Intended Use	14
4	Scope.....	15
4.1	Intended Use	15
4.2	Contraindication	15
5	Scope of delivery	15
6	Application	16
6.1	Functional check in the sterile operation area	16
6.2	Operating	16
7	Assembly and disassembly.....	17
7.1	Assembly of blade	17
7.2	Cleaning	18
7.3	Blade assembly	18
8	Reprocessing	19
8.1	Mechanical Reprocessing	19
8.2	Manual reprocessing	19
9	Storage	20
10	Warranty	20
11	Manufacturer and Service address	20

1 Directions for user manual

! In case of disregarding the user manual, there is the possibility of serious or fatal injury to the patient or operator!

Make sure that the following directions are understood and complied with!

- Every user is to read the manual completely and follow all instructions.
- Special attention must be paid to all information about precautions, warnings and dangers.
- The user manual must be readily accessible to the operator at any time.

2 Application and Safety

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

Wrong application and disregard of cautionary instructions can cause serious incidents. Therefore please inform yourself with the assistance of this user manual, as well as the instructions of the electrosurgical unit, about all functions and the safe handling of this instrument. Improper handling or maintenance, as well as inappropriate use can lead to premature wear and tear of the instrument.

- !** All instruments must be completely cleaned, disinfected and sterilized, before being put into use and before every other further application.
- !** It is very important to inspect each surgical instrument for visible damage and wear, e.g. cracks, fractures or fissures in the insulation before every use. In particular areas such as blades, points, locks, ratchets, latches and any moving parts, insulation and ceramic parts are to be carefully examined.
- !** Do not use damaged instruments.
- !** Do not use in the presence of flammable or explosive substances.
- !** The instrument should not be placed on the patient.
- !** Bear in mind the limitations of this treatment method! Do not attempt to seal vessels which are wider than the maximum diameter of 7 mm!
- !** Use the instrument solely according to its intended purpose.
- !** Coagulation of fallopian tubes with the BiSeal® instrument is **not** a safe method of sterilization!
- !** A functional check must be done before implementation.
- !** The blade carrier must be exchanged after each application. Never use the same blade carrier twice!
- !** When submerging the jaw parts into electrically conducting liquids, no safe vessel sealing is possible, due to the current being conducted by the liquid (shunt circuit)!

!	Even after ending the coagulation process, the outer surfaces of the jaw parts can be very hot and cause unwanted burns. Any touching of tissue immediately after the end of the sealing process is to be avoided!
!	Do not activate the blade carrier before the sealing of the grasped vessel is completed – danger of bleeding! Only after the electrosurgery unit has confirmed the proper vessel sealing process, the cut can be made.
!	Between the branches, in the blade duct and in the joints, rest substances might get deposited, which can cause a deterioration of function. (Electric shunt circuit, blade carrier gets stuck). In such cases the instrument must be cleaned during applications by use of delivered cleaning tools, a brush or a swab, as well as sterilized water.
!	During vessel sealing the blade must be retracted into the jaw parts, otherwise it can cause a short circuit and no sealing function can be effected.
!	Start with a low power output and increase the output according to requirement. If the sealing function is poor or not at all effective, look for cause. Do not simply turn the output power up to maximum.
!	The connection cable should be aligned in such a way that it does not touch the patient, nor has any contact with other cables. Protect the cable from mechanical damage (running over, squeezing, bending).
!	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
!	Avoid the simultaneous use of HF surgery and electrically grounded electrodes and sensors of diagnostic devices!
!	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
!	Observe manufacturer's instructions for electrical load capacity!


Additional directions for laparoscopic applications

!	Only do vessel sealing when the jaw parts of the instrument are visible for the operator. Unintentional activation or motion of instrument can lead to injury of patient!
!	The electrodes inside the jaw parts carrying current must have a safety distance of at least 2 cm to metallic end of trocar, when activated.
!	Insert and withdraw the instrument into/from the trocar with closed jaw parts only!

3 Intended Use

The vessel sealing instrument BiSeal® is intended to be used exclusively with units which are equipped with the BiSeal®-function from HEBU medical GmbH. Attention must be paid to the directions given in the manual instructions as delivered with the unit. The BiSeal® vessel sealing instrument is used for minimally invasive surgery, in particular for laparoscopy. The instrument can be inserted through the 10 mm trocar. Its purpose is to permanently seal veins, arteries and bundles of tissue, and thereafter implement the separation by cutting, without having to change instruments. Hereby vessels with a diameter of up to 7 mm can be permanently sealed. A cable connects the instrument to the HF-generator's bipolar sockets.

The vessel sealing instrument BiSeal® is available in working lengths of 200 mm (HF 9581-20) and 370 mm (HF 9581-37).

Due to interaction of mechanical pressure (p) and high-frequency current (I) a permanently and reliably sealed area of veins, arteries or tissue bundles is created. Hereby it is of significant importance that the pressure is being kept constant during the entire application period.

! Instruments for coagulation may only be used by persons who received special training or instructions in the use of such. Vessels with a diameter of up to 7 mm can be permanently sealed. For large vessels it is recommended to make a double coagulation, one next to the other. After sealing, it can be cut in the middle.

4 Scope

The use of HF bipolar clamps is used in HF surgery for coagulation and tissue cutting. In addition, for the permanent closure of veins, arteries and tissue bundles with subsequent separation.

4.1 Intended Use

HF coagulation forceps, laparoscopy; HF Clamps: A bipolar product covered with electrical insulation down to the working end. It is used in open surgery to coagulate vascular or tissue bundles for hemostasis. It is connected either to an electrosurgical diathermy unit or to an electrocautery generator that delivers the bipolar electrical current. It is typically designed as a scissors-type instrument with a long, curved, pincer-like mouth. It is made of metal and synthetic material in different sizes and shapes. It is a reusable product.

4.2 Contraindication

The RF bipolar clamps should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.

5 Scope of delivery

The instrument must be checked for correct contents as well as for being undamaged immediately upon delivery. In case of damages occurred during transport, these are to be reported right away.

HF 9581-20 or HF 9581-37	Vessel sealing instrument BiSeal® 200 mm or 370 mm
HB 2310	Needle holder – for save removal of the blades
HF 9581-01	Disposable replacement blades
-	Manual Instruction

6 Application

6.1 Functional check in the sterile operation area

The mechanical function must be checked before begin of operation, and by laparoscopic applications before inserting into the trocar: Opening, closing, locking, unlocking, moving the blade.

The electric function of the instrument should be checked before each application, e.g. on a sterile swab soaked with NaCl-solution.

6.2 Operating

After correct assembly and connection to one of the bipolar sockets on the HF-generator, set to BiSeal® mode. Then grasp the tissue which is to undergo treatment. Afterwards close the handle. The sealing current is activated by use of the foot switch.

The system ends the process by an automatic stop-function (automatic stop of flow of the high-frequency current, audible confirmation).


Following this, the sealed tissue can be cut with the integrated blade.

! Before beginning the sealing process, make sure that no ureter, nerve cords, intestinal loops or any other structures are grasped unintentionally!


! For large vessels it is recommended to make a double coagulation, one next to the other. After sealing, it can be cut in the middle.

7 Assembly and disassembly

7.1 Assembly of blade

1. Pull down safety latch and push halfway in


2. Remove blade with needle holder and dispose according to safety standards.

Attention!


Never remove blade with fingers – danger of injury!!!

7.2 Cleaning

3. Clean upper and lower parts of blade duct

7.3 Blade assembly

4. Insert new blade until it is caught in. Only then remove protective cap.

8 Reprocessing

8.1 Mechanical Reprocessing

We guarantee a minimum of 50 reprocessing cycles for the surgical instrument BiSeal®. The blade carriers are devised as a non-sterile single-use product. The warranty expires after using them once.

The durability of the instrument depends on its functions and on careful handling. Instruments for electro surgery are naturally subject to a higher wear and tear, depending on the way and duration of usage.

Preparation:

Remove coarse dirt from the instruments immediately after application. Do not use fixation agents or hot water (> 40 ° C), because this may lead to the fixation of residues and may impair the cleaning results.

!!! Important: also clean inside at thread of jaw parts!!!

Do not use any kind of acidly cleaning detergents. Substances with a neutral pH-value are most suitable.

Disinfection:

Mechanical, thermic disinfection should be carried out under consideration of the national guidelines with regard to the A0 value (see ISO 15883).

Drying:

Dry the outside of the instrument using the drying cycle of the cleaning/disinfection unit. Additionally, if necessary, it can be dried manually using a lint-free cloth. Dry the cavities of instruments with sterile air.

8.2 Manual reprocessing

Cleaning:

!!! Instruments must be disassembled before cleaning!!!

Prepare a cleaning bath according to manufacturer instructions.

- Rinse the instruments with cold tap water (<40 ° C) until all visible contaminants are removed. Imbedded dirt is removed with a soft brush.
- Submerge the instruments completely in the prepared bath. Induction period according to manufacturer's instruction.
- Clean the submerged instruments manually with a soft brush. All surfaces must be brushed several times.
- Rinse the products very thoroughly with running tap water for complete removal of detergent.

Disinfection:

- Prepare a disinfectant bath according to information provided by the disinfectant's manufacturer.
- Place the instruments into the disinfectant and observe the prescribed contact time.
- Rinse the products very thoroughly with deionized water to completely remove the disinfectant.

Drying:

The manual drying is effected by a lint-free cloth and, in particular for drying cavities and canals, by sterile compressed air.

Functional check:

Visual inspection for cleanliness, if applicable, assembly and functional test according to the manual. If necessary, repeat the treatment process until the instrument is visually clean. The original shine can be restored by polishing with a soft cloth.

Steam sterilization:

Sterilization of instruments by fractional pre-vacuum procedure (according to ISO 13060 / ISO 17665) taking into account the national guidelines of the respective country.

- 3 pre-vacuum phases with at least 60 mbar
- Heating to a sterilization temperature of at least 132 ° C, max. 137 ° C
- Shortest time: 3 min.
- Drying time: at least 10 min.

9 Storage

Store sterilized instruments in a clean, dry and dust-free environment at moderate temperatures of 5°C to 40°C.

10 Warranty

We provide a guarantee for all HEBUmedical products in accordance with statutory and country-specific legislation (proof of purchase by invoice or delivery note). The warranty period is two years.

Damage caused by incorrect handling or normal wear and tear is not covered by the warranty.

In case of warranty claims, please enclose a copy of the invoice or delivery note with the device.

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HEBUmedical GmbH