HF-Neutralelektroden, wiederverwendbar / HF return plates, reusable

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DEUTSCH

Operating Manual

ENGLISH

Inhalt

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Contents

- 1 Symbols 12
- 2 Introduction 12
- 3 Scope 12
 - 3.1 Intended Use 12
 - 3.2 Contraindications..... 12
- 4 Warning 13
- 5 Handling 13
- 6 Commissioning the unit..... 13
 - 6.1 Neutral electrode settings..... 13
 - 6.2 Application information 14
 - 6.3 Areas for attachment of the neutral electrode 14
 - 6.4 Connection to HF unit in monopolar mode..... 15
- 7 Hazard warnings 16
- 8 Protective measures and warnings..... 17
- 9 Instrument preparation 17
- 10 Maintenance, inspection 17
- 11 Packaging 17
- 12 Storage..... 17
- 13 Warranty / Repair 18
- 14 Manufacturer and servicing address..... 18

1 Symbols

| Symbol | Definition |
|---|--------------|
|  | CE marking |
|  | Attention |
|  | Manufacturer |

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

The use of HF neutral electrodes takes place in HF surgery.

3.1 Intended Use

HF Neutral electrodes: A non-active conductor attached to a patient and connected to a high frequency electrosurgical generator to establish a circuit for the return of electrical current to the generator after the transmission of the current for performing electrosurgical procedures on the patient. It is attached to the patient where the largest surface can be covered as close to the surgical field as possible. It is a reusable article.

3.2 Contraindications

The HF neutral electrodes must not be used on damaged and injured skin.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

4 Warning

| | |
|---|--|
| ! | Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use. |
| ! | A visual inspection of the instrument for damage and contamination must be carried out before each use! |
| ! | Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation! |
| ! | A visual inspection of the instrument for damage and contamination must be carried out before each use! |
| ! | The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return. |
| ! | Please observe the additional information enclosed with the products. |
| ! | Remove all protective sleeves and films prior to first using or preparation for use. |
| ! | The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user. |
| ! | Avoid improper throwing or dropping of instruments |
| ! | To avoid any damage, HF neutral electrodes with a damaged surface must be removed immediately! |
| ! | If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse. |

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Commissioning the unit

6.1 Neutral electrode settings

In the monopolar operating mode, a neutral electrode is required.

After attaching the neutral electrode to the patient, the HF unit recognizes the neutral electrode automatically (cf. illustration on the left and chapter **Fehler! Verweisquelle konnte nicht gefunden werden. Fehler! Verweisquelle konnte nicht gefunden werden.**). The HBS x-touch 200i recognizes neutral electrodes with divided and undivided contact surfaces.

If the displayed neutral electrode type does not agree with the used neutral electrode type, this is indicative of poor adhesion of the neutral electrode on the patient, a defective neutral electrode or a defective cable!

Please note the following information when attaching the neutral electrode

6.2 Application information


- The neutral electrode must be attached in such a way that its entire surface adheres to the patient's body. A rubber band provided for this purpose can be used to ensure better adhesion.
- When attaching the neutral electrode, ensure that the long side is pointing towards the operation field.
- After each change in the patient's position, check the correct adhesion of the electrode and the connecting cable!
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- The electrode must not be wet or coiled round.
- Fluids between the patient's skin and the neutral electrode must be avoided.
- Do not apply any additional electrically conductive gel to the neutral electrode.
- When performing surgery on small children, special electrodes suitable for children or babies must be used.
- When removing the neutral electrode, never pull on the cable!
- Never under any circumstances carry out repairs on the neutral electrode!

Please note that rubber electrodes lose their conductive capability if active substances are eroded due to cleaning of the material. This type of electrode considerably increases the risk of burns. Consequently ensure that not only the unit but also the reusable neutral electrode is included in regular maintenance inspections.


6.3 Areas for attachment of the neutral electrode

- Attach the neutral electrode to clean, dry skin.
- The electrode must be attached close to the operation area, but with a minimum direct distance from the operating field of 20 cm.
- Attach over tissue with a good supply of blood without dips or curves in the skin, for instance on the upper arm or thigh. (cf. picture below)
- Do not attach the neutral electrode to projecting surfaces.
- Do not attach above bones, scars, cuts or scratches.
- The electrode must not be attached to surfaces with a high amount of hair growth. If necessary, shave the relevant area before attaching.
- Do not attach above particularly adipose areas such as the abdomen or buttocks.
- The neutral electrode must not be attached above implants.


- Electrodes of monitoring systems must be placed at a distance of at least 20 cm from the operation area and from the position of the attached neutral electrode
- Correct attachment of the neutral electrode in adult patients

- Correct attachment of the neutral electrode in children

- Attachment in patients with pacemakers
Note! Consult the treating cardiologist!

RIGHT

WRONG!

6.4 Connection to HF unit in monopolar mode

First fix the neutral electrode ("return plate") to the patient, preferably on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. If a reusable neutral electrode (return plate) is used, do not apply conductive gel. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run in a transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the "active" HF cable to the HF unit.

7 Hazard warnings

The products may only be used in combination with HEBU accessories and HEBU equipment and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG). Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

For cleaning the active electrode, you must deactivate the "Auto Start" mode (only in bipolar applications) or disconnect the instrument from the generator. Directly after activation the tips could be hot. This may lead to burns.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

HEBUmedical GmbH assumes no warranty for the use of HEBU neutral electrodes in conjunction with HF devices and accessories from other manufacturers.

8 Protective measures and warnings

Use the greatest care when working with the neutral conductor!

Prevent the penetration of fluids between the patient's body and the neutral electrode.

A marked drop in power at a normal setting can be an indication of incorrect application of the neutral electrode.

9 Instrument preparation

Mechanical precleaning of neutral electrodes can be carried out). If a Washing-Desinfection unit is used for cleaning: max 95 °C.

The neutral electrode is **sterilized in the autoclave at 121°C / 20 min minimum dwell time**. When using the fractionated vacuum process, **sterilization is at least 134°C (USA 132°C) with a minimum holding time of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

Do not sterilize in hot air!

Cleaning should also be carried out immediately after the end of use in order to prevent residues from drying on.

10 Maintenance, inspection

After cooling to room temperature, the products must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

11 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container for sterilized packaging.

12 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

13 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

14 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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