HF-Bipolare Pinzetten

HF Bipolar Forceps / Pince bipolaire HF

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Seite 3

Operating Manual

Page 9

Mode d'emploi

Page 15

DEUTSCH

ENGLISH

FRANÇAIS

Inhalt

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Contents

1	Used Symbols	10
2	Scope	10
2.1	Intended Use	10
2.2	Contraindication	10
3	Storage	11
4	Cleaning and sterilization	11
5	Instructions for Use	11
5.1	Connection to HF unit in bipolar mode	11
6	Warnings	12
7	Safety Remarks	13
8	Servicing and manufacturer address	13

1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

2 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

The use of HF bipolar forceps takes place in HF surgery for the coagulation of tissue.

2.1 Intended Use

HF bipolar forceps, flex Tip: A bipolar product covered with electrical insulation down to the working end. It is used in open surgery to coagulate vascular or tissue bundles for hemostasis. It is connected either to an electrosurgical diathermy unit or to an electrocautery generator that delivers the bipolar electrical current. It is typically embodied as a long, tweezer-like instrument with tapered, flat jaw members that are closed by the surgeon by applying pressure to the connected legs. It is made of metal and synthetic material in different sizes and shapes. They are reusable products.

They are available with interchangeable tips. The flex-TIP system is available for replacing the forceps tips. The forceps tips are disposable products.

Maximum electrical capacity: 0,5kV_p.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore, carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

2.2 Contraindication

The HF bipolar forceps should not be used in conjunction with the central nervous or circulatory system.

3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

For cleaning, the cable, adapter, electrodes, electrode handles and shaft extensions must be disconnected from each other.

Thoroughly rinse the cable, connector, handle or electrode with deionized water.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

5.1 Connection to HF unit in bipolar mode

In bipolar HF surgery the coagulation is restricted to the area between the two electrodes of the bipolar instrument. Consequently, a neutral electrode is necessary. Otherwise, follow the recommendations regarding patient preparation in monopolar operating mode and the corresponding instructions for the HF unit. When using the bipolar "Auto start" mode, the cable used should not be longer than 4.5 m. Alternatively, use a twin cable with max. length of 5m.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

For cleaning the active electrode, you must deactivate the "Auto Start" mode (only in bipolar applications) or disconnect the instrument from the generator. Directly after activation the tips could be hot. This may lead to burns.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety Remarks

	Check the isolation on the electrode and the shaft extensions for damage before each use.
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	Do not rebend needles, loops or knives. Spatulas can be bent.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.

8 Servicing and manufacturer address

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HEBUmedical GmbH