# HF-Monopolare Klemmen

## HF monopolar clamps

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Seite 3

**Operating Manual**

Page 9

DEUTSCH

ENGLISH

# Inhalt

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# Contents

<b>1</b>	<b>Used Symbols</b> .....	<b>10</b>
<b>2</b>	<b>Scope</b> .....	<b>10</b>
2.1	Intended Use .....	10
2.2	Contraindication .....	10
<b>3</b>	<b>Storage</b> .....	<b>11</b>
<b>4</b>	<b>Cleaning and sterilization</b> .....	<b>11</b>
<b>5</b>	<b>Instructions for Use</b> .....	<b>12</b>
5.1	Connection to HF unit in monopolar mode.....	12
<b>6</b>	<b>Warnings</b> .....	<b>13</b>
<b>7</b>	<b>Safety Remarks</b> .....	<b>14</b>
<b>8</b>	<b>Servicing and manufacturer address</b> .....	<b>14</b>

## 1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

## 2 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

The use of HF bipolar clamps is used in HF surgery for coagulation and tissue cutting. In addition, for the permanent closure of veins, arteries and tissue bundles with subsequent separation.

### 2.1 Intended Use

HF coagulation forceps, laparoscopy; HF Clamps: A bipolar product covered with electrical insulation down to the working end. It is used in open surgery to coagulate vascular or tissue bundles for hemostasis. It is connected either to an electrosurgical diathermy unit or to an electrocautery generator that delivers the bipolar electrical current. It is typically designed as a scissors-type instrument with a long, curved, pincer-like mouth. It is made of metal and synthetic material in different sizes and shapes. It is a reusable product.

Maximum electrical capacity: 4,3kV<sub>p</sub>

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

### 2.2 Contraindication

The HF bipolar clamps should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.

### 3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

### 4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

The disposable monopolar electrodes and electrode handles are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Do not clean HF cables, adapters or electrode handles in an ultrasonic bath.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

Cleaning should also be carried out immediately after the end of use in order to prevent residues from drying on.

The instrument must be sufficiently dry after sterilization to avoid electric breakdown.

If a Washing-Disinfection unit is used for cleaning: max 95 °C. The HF cable or handle also may be immersed in a combined cleaning and disinfection solution (e.g. neodisher® MediClean forte, Dr. Weigert Co.). In such a case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Please observe the chemical manufacturer's instructions. Then rinse cable, connectors, and handle or electrode thoroughly with demineralized water.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

## 5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

### 5.1 Connection to HF unit in monopolar mode

First fix the neutral electrode ("return plate") to the patient, preferably on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. If a reusable neutral electrode (return plate) is used, do not apply conductive gel. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run in a transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the "active" HF cable to the HF unit. In order to avoid interference with camera system, do not run HF cables directly in line with camera cables. Never place HF cables directly on the patient's skin. Never lay HF cables in loops. Next, connect the instrument or handle with coagulation electrode, and switch on the HF unit. The electrode with the proper shaft diameter (4 mm or 2.4 mm) must be firmly attached to the handle tip. For handles with finger switches:

Yellow switch = actuation of cutting current

Blue switch = actuation of tissue coagulation

Always refer to the instruction manual of the HF unit and to the general information it provides relating to electrosurgical procedures.

## 6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

**Endogenous burns** are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

**Exogenous burns** are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

**Cardiac pacemakers** may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

## 7 Safety Remarks

	Check the isolation on the electrode and the shaft extensions for damage before each use.
	Do not rebend needles, loops or knives. Spatulas can be bent.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	Observe manufacturer's instructions for electrical load capacity!

## 8 Servicing and manufacturer address

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**Our vast range of products covers over 10 000 various instruments.  
They can be found in Internet or request our catalog.**

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HEBUmedical GmbH