HF- Bipolare Elektroden

HF bipolar electrodes

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Seite 3

Operating Manual

Page 9

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Used Symbols	10
2	Scope	10
2.1	Intended Use	10
2.2	Contraindication	10
3	Storage	11
4	Cleaning and sterilization	11
5	Instructions for Use	11
5.1	Connection to HF unit in bipolar mode	11
6	Warnings	12
7	Safety Remarks	13
8	Servicing and manufacturer address	13

1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer

2 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

The use of HF bipolar electrodes is used in HF surgery for coagulation, volume reduction and soft tissue cutting.

2.1 Intended Use

HF coagulation electrodes, ENT; HF Bipolar coagulation electrode, Laparoscopy; HF Bipolar electrodes, Arthroscopy; HF electrodes, special: A bipolar electrical conductor (blade or tip) connected to a mating electrical handpiece for establishing an active electrical connection between the output terminals of an electrosurgical generator and a patient. The delivered energy is used to perform lesions (cutting and / or cauterizing) at a desired location on the patient. It is commonly used in various medical disciplines. It is a reusable product.

Maximum electrical capacity: 1,0kV_p.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore, carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

2.2 Contraindication

The HF bipolar electrodes should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application.

3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

4 Cleaning and sterilization

Products are supplied in an unsterile state and must be cleaned and sterilised before use.

The disposable monopolar electrodes and electrode handles are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Do not clean HF cables, adapters or electrode handles in an ultrasonic bath.

Do not sterilise HF cables, adapters or electrode handles and electrodes in hot air sterilisers.

Cleaning should also be carried out immediately after the end of use in order to prevent residues from drying on.

The instrument must be sufficiently dry after sterilization to avoid electric breakdown.

If a Washing-Disinfection unit is used for cleaning: max 95 °C. The HF cable or handle also may be immersed in a combined cleaning and disinfection solution (e.g. neodisher® MediClean forte, Dr. Weigert Co.). In such a case, the service life may then be reduced due to oxidation of the contacts and crystallisation in the plugs. Please observe the chemical manufacturer's instructions. Then rinse cable, connectors, and handle or electrode thoroughly with demineralized water.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

5.1 Connection to HF unit in bipolar mode

In bipolar HF surgery the coagulation is restricted to the area between the two electrodes of the bipolar instrument. Consequently, a neutral electrode is necessary. Otherwise, follow the recommendations regarding patient preparation in monopolar operating mode and the corresponding instructions for the HF unit. When using the bipolar "Auto start" mode, the cable used should not be longer than 4.5 m. Alternatively, use a twin cable with max. length of 5m.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

For cleaning the active electrode, you must deactivate the "Auto Start" mode (only in bipolar applications) or disconnect the instrument from the generator. Directly after activation the tips could be hot. This may lead to burns.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety Remarks

	Check the isolation on the electrode and the shaft extensions for damage before each use.
	Do not rebend needles, loops or knives. Spatulas can be bent.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	Observe manufacturer's instructions for electrical load capacity!

8 Servicing and manufacturer address

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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HEBUmedical GmbH