HF- Monopolare Elektroden, Einweg

HF monopolar electrodes, disposable

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Seite 3

Operating Manual

Page 9

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Used Symbols	10
2	Scope	10
2.1	Intended Use	10
2.2	Contraindication	11
3	Storage	11
4	Cleaning and sterilization	11
5	Instructions for Use	11
5.1	Connection to HF unit in monopolar mode.....	11
6	Warnings	12
7	Safety Remarks	13
8	Servicing and manufacturer address	13

1 Used Symbols

Symbol	Definition
	CE-marking
	Attention
	Keep away from sunlight
	Store at dry place
	No intervention on patients with cardiac pacemakers.
	Manufacturer
	Do not reuse
	Sterile

2 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

The use of monopolar disposable electrodes is used in HF surgery for coagulation and tissue cutting.

2.1 Intended Use

HF disposable monopolar electrodes: A monopolar electrical conductor (blade or tip) connected to a mating electrical handpiece for establishing an active electrical connection between the output terminals of an electrosurgical generator and a patient. The delivered energy is used to perform lesions (cutting and / or cauterizing) at a desired location on the patient. It is commonly used in various medical disciplines. The electrode can be monopolar or bipolar and is a disposable product. Maximum electrical capacity: 4,3kV_p.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore, carry out a thorough visual examination and a function test prior to each application of these products. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

2.2 Contraindication

The HF monopolar electrodes should not be used in conjunction with the central nervous or circulatory system. Tube ligation is not a safe application. Patients with cardiac pacemakers and in-vivo cardiac defibrillators have limitations in their use. RF Monopolar Electrodes should not be used in the immediate vicinity of bone, tooth root, and metal.

3 Storage

Temperature: -20 °C - +50°C

Relative humidity: 0 -75%, not condensing

Atmospheric pressure: 500 -1600 mbar

Products must be protected from direct sunlight in storage. The storage of the products in their original packaging until their first use is recommended.

4 Cleaning and sterilization

STERILE

The disposable monopolar electrodes and electrode handles are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

This product is a disposable product and must **NOT** be reused!

5 Instructions for Use

Thoroughly examine all products for visible irregularities. Do not use cables with a brittle or defective insulation. In case of doubt check operation by means of a suitable testing device.

5.1 Connection to HF unit in monopolar mode

First fix the neutral electrode ("return plate") to the patient, preferably on the upper arm or on the thigh. Ensure that the skin is free from hair and grease. If a reusable neutral electrode (return plate) is used, do not apply conductive gel. The current paths inside the patient's body should be short and must run diagonally. Current paths must never run in a transverse direction through the body or across the thorax. The patient must be insulated against all electrically conductive parts. Ground the operating table. Place the patient on a dry electrically insulating layer. Prevent skin areas touching each other. Place dry gauze in between. Then connect the neutral electrode, the foot switch and the "active" HF cable to the HF unit. In order to avoid interference with camera system, do not run HF cables directly in line with camera cables. Never place HF cables directly on the patient's skin. Never lay HF cables in loops. Next, connect the instrument or handle with coagulation electrode, and switch on the HF unit. The electrode with the proper shaft diameter (4 mm or 2.4 mm) must be firmly attached to the handle tip. For handles with finger switches:

Yellow switch = actuation of cutting current

Blue switch = actuation of tissue coagulation

Always refer to the instruction manual of the HF unit and to the general information it provides relating to electrosurgical procedures.

6 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

7 Safety Remarks

	Check the isolation on the electrode and the shaft extensions for damage before each use.
	Do not rebend needles, loops or knives. Spatulas can be bent.
	If damaged, do not use these products!
	Do not use in the presence of combustible or explosive materials!
	Never lay instruments on the patient or in his/her direct vicinity!
	HEBUmedical expressly warns against modifying the product. Any modification exempts HEBUmedical in from any and all liability. .
	Protect these products from any form of mechanical damage! Do not throw! Do not: use force!
	Defective HF cables or electrode handles may not be repaired. Always replace a defective part with a new one.
	Observe manufacturer's instructions for electrical load capacity!

8 Servicing and manufacturer address

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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HEBUmedical GmbH