HB 0181-06

HB 0183-01

Stimmgabel / Tuning fork

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Seite 3

Operating Manual

Page 7

DEUTSCH

ENGLISH

HB 0181-06
HB 0183-01

Inhalt

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HB 0181-06
HB 0183-01

Contents

1	INTENDED USE	8
2	APPLICATION NOTES	8
3	CLEANING.....	8
4	RYDEL-SEIFFER-TUNING FORK	8
4.1	HB 0183-01 Tuning fork with a C 64 damper and a c 128 base.....	8
4.2	HB 0181-06 Tuning fork with damper and base, adjustable from c 128 to h 240.....	9
5	MANUFACTURER AND SERVICING ADDRESS	10

1 Intended use

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

A tuning fork is a u-shaped instrument, usually made of stainless steel, with a handle at the base of the 'u'. The vertical prongs of the 'u' are cut to a particular length to produce sound of a specific wavelength when the 'u' is struck against a hard object, usually rubber. It is typically used during ear/nose/throat (ent) examinations to test hearing acuity. This is a reusable device.

2 Application Notes

- If possible, the tuning fork should be struck against the ball of the thumb or against a piece of wood or plastic. Never strike it against hard objects (e.g. metal, steel or stone). It is important that the striking is short and bouncy i.e., the tuning fork is touched only briefly. The period of vibration of the sound achieves a maximum period of time by this short striking.
- Now the required examinations can be performed optimally.
- If the tuning fork has fixed weights, they should not be removed. The frequency changes by moving the weights and an exact adjustment can only carry out by the manufacturer.

3 Cleaning

After use, the instrument should not be cleaned in sterilizers. To clean tuning forks, they should be rubbed off with a dry or lightly oiled cloth.

4 Rydel-Seiffer-Tuning fork

4.1 HB 0183-01 Tuning fork with a C 64 damper and a c 128 base

The Rydel-Seiffer Tuning fork produces two different frequencies: note C = 64 Hz with the dampers (as marked on the dampers) and note c = 128 Hz without the dampers (as marked on the tuning-fork itself). The dampers must be fitted with the lower edge flush with the gradation on the prongs. The imprint „C 64,, on the dampers, the imprint „c128,, on the front of the tuning-fork and the gradation must all be facing the user when the damper is fitted. The knurled-head screws are tightened to fix the dampers and slackened to remove them. They should be tightened firmly when the dampers are fixed.

On examination (sensitivity test, etc.):

After striking with the black the tuning fork will be placed on the points of the patients, which will be examined.

On the back side of the tuning fork scales are fixed on the dampers. By observation of these scales, the intensity of vibration of the tuning fork can begin each with the same intensity of vibration.

This is done as follows:

The tuning fork will struck against the ball of the thumb with the dampers (the sound C=64z). Initially, by the strong vibration of the prong of the tuning fork both triangles are blurry to see.

After a few seconds, depending on the striking strength, the triangles start to be clearly visible at the baseline. In addition, they grow slowly, that means with slowly decreasing amount of the intensity of vibration until the first, second, third, etc. cross of the scale. So, the examination can be started each at a particular cross of the scale.

This can be read off either the white or the black triangle, depending on the amount of light, available or which scale is more easily readable.

4.2 HB 0181-06 Tuning fork with damper and base, adjustable from c 128 to h 240

The Lucae-tuning fork produces 7 different frequencies, from c 128 to h 240. The exact frequencies are listed in the following table:


	c	d	e	f	g	a	h
Frequence in Hz	128	144	160	170,66	192	212,33	240

The change of the frequency will carry out by moving the weights at the prongs. These weights can be adjusted manually.

When the sound is changing, the spring in the weight have to lock in the groove of the tuning fork. Only therefore an exact frequency can be reached.

Both weights have to be always on the same level. Are the weights on different levels, it's not possible to produce an intermediate sound.

The stamped letter directly under the weight identifies the produced sound.

5 Manufacturer and servicing address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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