**HB 1269-00**

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
# HB 1269-00

## LED High performance cold light source

### Content

- 1 Device overview..... 16
- 2 Important notes..... 18
  - 2.1 Symbols ..... 18
  - 2.2 Designates use ..... 18
  - 2.3 Safety notes ..... 18
- 3 Scope ..... 20
  - 3.1 Inteded Use..... 20
  - 3.2 Contraindications ..... 20
- 4 Operation..... 20
  - 4.1 Commissioning..... 20
  - 4.2 Brightness control ..... 20
  - 4.3 Display ..... 21
  - 4.4 Filters (not available for all models)..... 21
  - 4.5 Light guide adapter (not available for all models)..... 22
  - 4.6 Module (optional) ..... 23
  - 4.7 Menu ..... 23
  - 4.8 Footswitch..... 24
  - 4.9 Remote control via USB..... 24
  - 4.10 Technical data..... 24
  - 4.11 Disposal ..... 25
- 5 Warranty..... 25
- 6 Manufacturer and Service address ..... 25


## 1 Device overview

### FRONT VIEW LED3000


- 1 Main switch (ON/OFF)
- 2 Rotary knob with pushbutton function
- 3 Display
- 4 Light guide connection
- 5 Filter holder
- 6 Handle

Design depending on device

**DEVICE REAR LED3000**

1) max. screw-in depth 7mm measured from outer housing surface  
 2) measured to vertex of front

- 7 ESD socket
- 8 Socket for footswitch
- 9 USB connection
- 10 Supply socket
- 11 cover for optional modules
- 12 light guide fixture
- 13 adapter fixture

## 2 Important notes

### 2.1 Symbols

SYMBOL	DESCRIPTION
	Warning of a danger area! Caution, please refer to documentation!
	Warning of dangerous electrical voltage!
	Warning of hot surfaces!
	ATTENTION: Dangerous optical radiation! Do not look into the light!

### 2.2 Designates use

This product has been developed and is intended for illumination purposes in the sector of stereo microscopy and endoscopy for industry and laboratory use. The applicable standards and guidelines are listed in the Chapter "Technical data".

This light source generates light that can be passed through a light guide to the microscope or object to be illuminated. Numerous light guides, such as e.g. flexible light guides, goosenecks or ring lights are available. Comprehensive accessories such as e.g. filters, focusing lenses and adapters are also available.

We would like to say thank you for choosing to purchase this high-quality product and wish you lots of success whilst using it.

### 2.3 Safety notes

	<p>The operating instructions must be read before using the device. It is recommended to store the operating instructions close to the device to ensure that it can be read immediately by the user in case of need. The user must inform a responsible person or the manufacturer immediately if there are any changes in the operating behaviour or to any safety-relevant device components. The operator is responsible for the correct use of the product and for the training of the user. The generally application legal national safety and accident prevention regulations must be complied with in addition to the supplementary user guidelines of the operator. The manufacturer is not liable for any damages due to incorrect operation or non-compliance with the danger instructions!</p>
	<p>Never open the device or its components, unless specifically indicated in this manual as you could otherwise risk physical injury or damage. If necessary, contact your retailer. Technical modifications to the device are not permitted under any circumstances and repairs should only be carried out by the manufacturer or authorised service points. Only original spare parts must be used. Only use original spare parts approved by the manufacturer for this product.</p>

	<p>Do not operate the device in the vicinity of water or in rooms with extreme humidity. Do not spill any liquids onto the device and do not push any foreign bodies into the product, you could be electrocuted or cause damage to the device. Do not position any flammable materials in the direct vicinity of the device. The device must be kept away from all sources of heat and is not approved for use in explosion hazard areas! The device may only be used in closed rooms.</p>
	<p>Never allow the device to operate unsupervised! The device may only be operated with the specified supply voltage. Only use as much brightness as required for working. A higher brightness setting can irritate the eyes, stresses the object and also reduces the service life of the device.</p>
	<p>Always keep the ventilation openings of the device free to avoid causing any overheating. Ensure that the free air circulation is not hindered by anything and that there is a sufficiently large gap between the ventilation openings and other objects. Keep the environment of the device clean (above all, dust and lint-free), so that any device fans cannot suck in any dirt.</p>
	<p>ATTENTION: dangerous optical radiation, do not look into the light, can injure your eyes</p>
	<p>The device emits light with a high intensity! Never look directly into the light, any light-emitting device openings or into the light guide end as you will risk injury to the eyes. The emitted light is absorbed by materials and converted into heat. Never cover the light emission aperture of the device or the light guide outlet end with the hand, any other parts of the body or any materials. To illuminate objects correctly, always ensure there is an appropriate gap between the object and the light guide outlet end, and an appropriate brightness level to avoid thermal damage or fire. Never leave the light guide end lying around freely.</p>
	<p>Before cleaning, disconnect the product from the mains and allow it to cool for at least 15 minutes. Use a dry cloth for cleaning or a cloth moistened with water or a 5% ethyl or isopropyl alcohol solution, but do not use combustible or flammable liquids under any circumstances. Allow the device to dry thoroughly before use. If any liquid penetrates the device accidentally, pull the plug immediately and bring the device to a service point. If dust has collected on the ventilation openings, it must be removed carefully. Excessively intensive cleaning or the use of solutions that are too concentrated can lead to scratches, discoloration or other damage to the surfaces of the device.</p>
	<p>Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.</p>

## 3 Scope

The use of light sources takes place in diagnostics.

### 3.1 Inteded Use

Cold light Source and Accessories: An electrically operated device (AC) used in general surgery or to examine a patient and designed to generate high intensity light, often called cold light. This light is usually conducted through a fiber optic cable that supplies a device with light. For example, headlights, microscopes and endoscopes are supplied with light during the treatment of a patient. There are special light sources for endoscopes for this purpose.

### 3.2 Contraindications

There are no known contraindications.

## 4 Operation

### 4.1 Commissioning

	<p>Unpack the light source and remove any protective foils. Some models are shipped with a protective cover at the light guide connection (4). Please remove it before turning on the light source and keep it for later use during transportation and storage.</p> <p>Pick up the light source with the handle (6) and position it on the required location. Insert the light guide to be used in the appropriate opening (4) on the light source and position the guide as required.</p>
	<p>Please ccheck the light guide for cleanness before inserting into the light source, otherwise you risk to burn the light guide and the light source. Due to the high light flux, glued fibers can be destroyed. Only use certified and tested fibers.</p>
 	<p>Please comply with the safety instructions while doing so!</p>
	<p>Connect the light source at the inlet (10) with the supplied power unit. With some power supplies the handle piece of the DC plug must be pulled backwards to fit correctly to the socket (10). Also pull handle piece when unplugging from (10). An appropriate ESD discharge cable can be connected to the ESD port (7) in ESD-sensitive areas. The light source can now be switched on by pressing the switch (1).</p>



### 4.2 Brightness control

Turn the rotary knob (2) to adjust the brightness in 100 steps. Please note that the light source is not switched off when the position "0" is selected.

Pressing the rotary knob (2) directly selects steps previously saved in the menu.

### 4.3 Display

The display (3) shows all information relevant to operation. It displays the current brightness level and also a graphic progress bar. Numerous symbols also provide information about the operating state. The meaning of these symbols is explained in the following table.

	Insert the light guide in the connection
	Prevents brightness control directly at the device
	Brightness control via remote
	Light source is connected to the PC via USB
	Standby (triggered by menu, footswitch or USB)
P1...P10	Displays brightness presets that have been called
	Progress bar, displays the brightness graphically

### 4.4 Filters (not available for all models)

Various filters can be installed as required in the filter holder (5). Simply pull the filter holder out of the device for this purpose. Only use clean filters approved by the manufacturer.

Please note that the filter can be hot during operation and immediately afterwards!

The following filters are available:

Item no.	Item description
HB1269-02	Plug-in filter for HB1269-00, halogen 3.500K
HB1269-03	Plug-in filter for HB1269-00, blue
HB1269-04	Plug-in filter for HB1269-00, green
HB1269-05	Plug-in filter for HB1269-00, yellow
HB1269-06	Plug-in filter for HB1269-00, red
HB1269-07	Plug-in filter for HB1269-00, daylight
HB1269-08	Plug-in filter for HB1269-00, grey

## 4.5 Light guide adapter (not available for all models)

Some models allow the use of adapters, which can be placed into the light guide connection (4) until it snaps. For hard fixing the adapter, tighten adapter fixture (13) with allen key.


For installing the light guide, loosen (12) and insert light guide into (4). While pressing smoothly the light guide into the connector (4), gently tighten (12) again. Check that screw (12) really is fixing the light guide and is not irradiating by the light.

Some light guide - adapter combinations allow a snap-in installation of the light guide without fixing the light guide with (12). In such case the screw (12) can be replaced by a cover, but store the screw (12) in a safe place for later use.

Please bear in mind that the light guide tip must be in close contact with the optics inside the light source for maximum light output. Tolerances can cause a small gap between light guide and optics when using snap-in light guide installation, which lead to losses in luminous flux.


Always Disconnect device from mains when installing an adapter. Please notice, that adapters can be hot, so let them cool down before changing.

## 4.6 Module (optional)

On the back at (11) special modules can be installed. E.g. an analog module for automated control of the light source via DIN sockets.

## 4.7 Menu

Menüpunkt	Wertebereich	Bemerkung
SETTINGS STARTUP LOWER LIMIT UPPER LIMIT PRESET 1-10 <MAIN MENU	0-100%, LAST 0-100% 0-100% 0-100%, SKIP	Brightness value after switch-on Smallest settable value Largest settable value SKIP Memory presets for brightness values, this memory preset is not used if SKIP is selected Back to the main menu
STROBE LEVEL PULSE RATE RUN/STOP <MAIN MENU	0-100%, LAST 0,1-5000, 0ms 0,2-5000, 0Hz ACTIVE, STOPPED	Brightness value for the strobe function Pulse duration Pulse frequency Start/Stop strobe function Back to the main menu
SYSTEM SHUTTER USB Mode  STEP SIZE FOOTSW.  FOOT DIM SPEED  PANEL LOCK FIRMW ARE RESET SETTINGS LCD CONTRAST LCD BRIGHTNESS <MAIN MENU	LIGHT ON, LIGHT OFF STANDARD LEICA 1-5% TOGGLE+DIM  TOGGLE PUSH ON PUSH OFF PRESETS DUAL KEY 1-10X  YES, NO  YES, NO, DONE 0-100% 0-100%	The shutter function can be used Photonic protocol Leica LAS protocol Step size for dimming Switch on/off with footswitch, hold footswitch down for dimming. Switch on/off with footswitch Light source activated while pressed Light source deactivated while pressed Press to scroll through Presets Dimming with 2 external key buttons Dimming speed for TOGGLE+DIM and DUAL KEY Prevents brightness control directly at the device Version display of the firmware Reset to factory settings Display contrast Display brightness Back to the main menu

The menu can be exited at any time by keeping the rotary knob pressed down.  
For new features please check our website or ask your local dealer.

## 4.8 Footswitch

In the standard setting, the light source can be shuttered, i.e. switched on/off by means of a footswitch connected to the relevant port (8). In addition, keeping the footswitch pressed down changes the brightness (dimming). Once the required brightness is reached, simply release the footswitch. Other functions can also be programmed depending on the menu setting.

## 4.9 Remote control via USB

The light source can be controlled from the PC via USB (9). You will need a commercial USB cable and a driver. All protocol information can be obtained on request from the manufacturer.

## 4.10 Technical data

Menüpunkt	LED 3000
Dimensions (WxDxH)	approx. 170x196x98 without protruding parts
Weight	approx. 3.8kg
Illuminants	High Power LEDs
Brightness	approx. 470lm with fiber Ø 5x1000mm approx. 640lm with fiber Ø 8x1000mm
Color temperature	approx. 5800K ±5%
Maximum fiber diameter	9mm
Service life L <sub>70</sub>	approx. 30,000h (70% of initial brightness)
Cooling system	Axial ventilator
Supply	12VDC, 5420mA
Power consumption	Max. 65 Watt
Operating conditions	10-40°C, max. 80% relative humidity, max. 2000m sea level
Storage and Transportation	-25 to +70°C, max. 80% relative humidity non. condensing
Light guide connection	Diameter 15mm
Connections at the back	USB: for control with PC Mono jack 2.5mm: Footswitch (tip: signal, back: ground) ESD socket DC socket
Display	Graphic display for operation display and comprehensive menu functions
Fixation	2xM5 thread in distance of 90mm
Approval (CE)	EMC directive 2004/108/EC
Classification EN62471:2008	Risk group 2
Mains wire	original or VDE/HAR 3x0,75mm <sup>2</sup> 10A/250V

## 4.11 Disposal

Do not dispose of the product in household waste, take it to the appropriate collection and disposal points. Comply with the national and local disposal guidelines and laws. Please help to keep our environment clean.

The development of this product is subject to continuous development. Although we try to keep this operating manual up to date, we reserve the right to change device data at any time without prior notice. The operating manual is produced subject to errors and printing errors.

## 5 Warranty

We provide a guarantee for all HEBUmedical products in accordance with statutory and country-specific legislation (proof of purchase by invoice or delivery note). The warranty period is two years.

Damage caused by incorrect handling or normal wear and tear is not covered by the warranty.

In case of warranty claims, please enclose a copy of the invoice or delivery note with the device.

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HEBUmedical GmbH