Standard-Instrumente / Standard Instruments

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Seite 4

Operating Manual

Page 22

DEUTSCH

ENGLISH

Inhalt

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Contents

1	Symbols	22
2	Introduction	22
3	Scope	22
3.1	Intended Use	24
3.2	Contraindications.....	29
4	Warning.....	30
5	Handling	30
6	Preparation.....	31
6.1	Reutilization restrictions	31
6.2	Information on instrument preparation	31
6.3	Preparation at the place of use	31
6.4	Ultrasound bath (optional)	32
6.5	Manual cleaning	32
6.6	Mechanical cleaning.....	32
6.7	Drying.....	34
7	Maintenance, inspection	35
7.1	Functional check	35
8	Sterilization.....	36
8.1	Packaging.....	36
9	Lifetime.....	36
10	Storage.....	36
11	Warranty / Repair	37
12	Manufacturer and servicing address.....	37

1 Symbols

Symbol	Definition
	CE marking
	Attention
	Validated Parameters
	Manufacturer

2 Introduction

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

3 Scope

The use of Standard Instruments is mainly in the following areas:

Gall Stone Probes: The instrument is used in surgical procedures of the gallbladder and bile duct. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Bile Duct Dilators: The instrument is used in surgical procedures of the bile duct. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Biopsy curettes, Uterine curettes; Sharp Hooks; Myoma Screw: The instrument is used for surgical procedures in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Eye Hooks; Strabismus Hooks; Eye Specula; Exenteration spoon; Chalazion spoon; Eye Needles: The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Lachrymal Dilators; Lachrymal Probes: The instrument is used for procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Septum and Periostic Elevators; Ring curettes; Elevator; Elevators: The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Lung spatulas: The instrument is used for procedures in thoracic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Curette; Currettes, titanium; Lupus currettes, -spoons, -scrapers; Bladder Retractors; Gall Stone Scoops; Uterine Currettes; Sharp spoon; Rectal director: The instrument is used in gastroenterological and urological surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Elevatorium: The instrument is used for various surgical procedures on the skull. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Hooks; Exploration-& Coagulation Hooks; Root retractors; excavators; excavator- curette; Tissue spatula Ligature Guides; Wire strut guide; Tissue Hooks; Micro Hooks; Spatula; Intestinal Needles; Muscular- & Skin suture; Reverdin Needles; Ligature Needles; post mortem needle; Guide needles; Ligature Conductors; Laser Instruments; Hooks; Extension Bows: The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Abdominal-& Intestinal Spatulas; Organ spatulas: The instrument is used in various surgical procedures in the abdominal area. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Tissue dissectors: The instrument is used for various procedures in general or plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Vascular Dilators; Endarterektomiestripper; Varicose vein cutlery; Venous Instruments, Venous Hooks, Venous Extractors: The instrument is used in various procedures in vascular surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Ear currettes; Paranasal Sinus Currettes; Hooks for posterior nares; Ear Hooks; paracentesis Needles; Foot plate instruments; Measuring probes; Ear Loops; Thyroid Gland retractors; Hooks; Micro-laryngeal Instruments; Probes; Ear Specula; Tongue spatulas; Comedo Extractors; Laryngeal Mirrors; Laryngeal mirror handles; Middle Ear Mirror; Cotton Applicators: The instrument is used during ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Areola markers; Breast Enlargement Sets: The instrument is used in plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Cement spatula; Table top cutters: The instrument is used in orthopedics. The application may only be carried out by suitably trained and qualified specialist personnel.

Uterine depressors; Uterine-Dilators: The instrument is used during examinations in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Prostatic Dilators; Dilatationsbougies: The instrument is used during examinations in urology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Plaster Spatula; Mallets; Screwdriver; Lead Hands, Filling instrument: The instrument is used in various procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

3.1 Intended Use

Curettes:

Curette; Curettes, titanium; Lupus curettes, -spoons, -scrapers: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon, secretion spoon. It is a reusable instrument intended for transient use.

Biopsy curettes, Uterine curettes: A heavy, hand-held, surgical instrument with an elliptically windowed end on a long shaft for scraping out the uterine lining, scraping tissue during a gynecological procedure, and removing waste from the uterus. The instrument is usually made of stainless steel and can be either rigid or flexible. The working end is prepared sharp or blunt in different versions. The general product group includes uterine, placental and fallopian tube curettes. It is a reusable instrument intended for transient use.

Uterine curettes: A heavy, hand-held surgical instrument with an elliptically windowed working end on a long shaft for scraping out the lining of the uterus, scraping out tissue during a gynecological procedure, and removing abortion residues from the uterus. The instrument is usually made of stainless steel and can be either rigid or flexible. The working end is sharp or blunt in various designs. The general product group includes uterine, placenta, fallopian tube curettes. It is a reusable instrument intended for transient use.

Ear curettes; Paranasal Sinus Curettes: A surgical hand instrument having an annular, fenestrated, or spoon-shaped working portion which may be either sharp or blunt to provide tissue by scraping or scraping during an ENT procedure. It is a curette for general use in ENT surgery, it is made of stainless steel and is available in various designs and shapes. It is a reusable instrument intended for transient use.

Ring curettes: A surgical hand tool with an annular windowed working part for use in orthopedic surgery. It is used to remove bone tissue. It is a reusable product.

Exenteration scoops; Chalazion Spoon: An ophthalmic surgical instrument with a windowed, spoon-shaped or ring-shaped end that is either sharp or blunt. It is used to scrape out eye tissue. It is usually made of stainless steel. It is a reusable instrument intended for transient use.

Exenteration spoon; Chalazion spoon: An ophthalmic surgical instrument with a windowed, spoon-like or ring-shaped working end that is either sharp or blunt. It is used to scrape out eye tissue. It is usually made of stainless steel. It is a reusable instrument intended for transient use.

Elevators:

Elevator; Elevators: A surgical hand instrument for lifting, positioning or levering of bone structures, other anatomical structures or surgical material during an orthopedic procedure. It is made of stainless steel and is available in a variety of sizes, styles and working ends. It is usually a sturdy instrument with a blunt working end that can be pointed, rounded, flat, curved or hooked. The handle at the proximal end is also available in different versions. Some models have special functions, e.g. for ankle surgery, or to treat a dislocation. It is a reusable instrument intended for transient use.

Septum and Periostic Elevators: A surgical hand instrument for lifting, positioning or rupturing anatomical structures or surgical material in the ENT area. It can be a single or double instrument with a blunt end (leaf or hook in different versions) and a handle in the middle. Size and design of intended use, e.g. a dull impeller elevator, a septum elevator, a stronger nose artery elevator, a special viable ear elevator. The instrument is usually made of stainless steel. It is a reusable instrument intended for transient use.

Elevatorium: A surgical hand instrument for levering out, lifting and positioning a bone plate during surgery on the skull. Also known as bone flaps Elevator there is typically made of one piece and is robust, long, rounded, with a blunt blade at the distal working end and a handle at the proximal working end. It is typically made of high quality stainless steel. It is a reusable instrument intended for transient use.

Tissue dissectors: A surgical hand instrument, usually in the form of spoon or rounded at the working end. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. The instrument is reusable.

Comedo Extractors: A hand instrument for levering out, lifting and positioning anatomical structures or surgical materials during ear, nose and throat (ENT) procedures. It may be a single- or double-ended instrument with blunt working end or hooks in various shapes without a cutting edge at the distal working end (s), and a handle positioned centrally or at the proximal end. Its shape and size depends on the intended use, for example: a blunt tympanic elevator, a rhinological or septal elevator, a strong nasal fracture elevator, or a specially shaped ears elevator. It is usually made of stainless steel. It is a reusable instrument.

Eye Specula: A surgical ophthalmic instrument for holding the edges and associated tissue of a cut in the eye during a surgical procedure. It usually has a self-retaining arm with curved hooks or leaves at the distal end that grip and hold the tissue. The proximal end is usually connected and has a lockable mechanism to hold the working end in an expanded position. It is made of stainless steel and is available in different sizes and designs. It is a reusable instrument intended for transient use.

Uterine depressors: A hand-held instrument for holding off the uterus in order to facilitate an examination of the surrounding organs / tissues. It is typically a double-ended product in various designs (e.g., curved, spoon-shaped, or fenestrated) with a slightly angled working end to better manipulate the organs. It is made of stainless steel and is a reusable instrument.

Hooks:

Tissue Hooks; Micro Hooks: A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks that apply tension to the skin. It is a reusable instrument intended for transient use.

Hooks; Exploration-& Coagulation Hooks; Root retractors; excavators; excavator- curette: A shaft-like surgical instrument, having a proximal end forming a handle which can be formed in various embodiments. The shaft at the distal end merges into one or more parallel hooks to hold tissue. It is a reusable instrument intended for transient use.

Thyroid Gland retractors; Hooks: A surgical instrument for stably holding the trachea during a tracheostomy and / or separating / pushing away the margins of a tracheostomy to set a tracheostomy tube. The instrument may be of various types, e.g. Single or double ended, with a wide retractor on one side and opposite a pointed hook or several parallel hooks, which can be bent at different angles. It is made of stainless steel. It is a reusable instrument intended for transient use.

Hooks for posterior nares: A shaft-like surgical instrument whose proximal end forms a handle that can vary in shape. The shaft tapers toward the distal end toward a hook or hooks used in plastic nose surgery. In the multi-hook variant, the hooks are arranged in parallel. It is a reusable instrument intended for transient use.

Ear Hooks; paracentesis Needles; Foot plate instruments: A fine surgical instrument that consists of a handle that tapers into a shaft and ends in a fine, 90 ° bent tip. It is used during ENT surgical procedures to repair structures in the middle ear, e.g. to edit the stapes. It is a reusable instrument intended for transient use.

Eye Hooks; Strabismus Hooks: An ophthalmic surgical instrument with a shaft-like handle in various versions that tapers towards the working end. The distal part of the instrument is blunt or pointed bent or angled to better manipulate the ocular tissue or to remove debris. It is a reusable instrument intended for transient use.

Laser Instruments; Hooks: An instrument used in conjunction with a suitable endoscope during endotherapy. It is used for mechanical work, e.g. grasp tissue or foreign objects, to bring medical fluid or width of an area / lumen. It works without electricity, including e.g. High-frequency, electromagnetic, ultrasonic or laser energy. It is a non-active product that is reusable and intended for transient use.

Sharp Hooks: A shaft-like surgical instrument whose distal end is in the form of a hook or more parallel hooks for fixing or pulling uterine fibroids during a gynecological procedure. It is typically a one-piece stainless steel instrument with a handle in various designs at the proximal end. It is a reusable instrument intended for transient use.

Dilators:

Uterine-Dilators: A round rod-like massive instrument that is inserted through the cervix and dilates the cervical canal. Commonly known as a uterine dilator, it is typically made of metal or plastic and has a tapered rounded distal tip. It can have one or two different diameters, one at each end of the instrument, and can be part of a set of size-increasing diameters which are inserted successively. The instrument allows insertion of a curette or other instrument into the uterus and / or expands the uterus by pushing against the walls of the uterus from the inside. The instrument is reusable.

Bile Duct Dilators: A surgical instrument consisting of a slender hollow or solid body of metal or plastic in a cylindrical shape in a variety of sizes for dilating the bile duct, e.g. at the union of the ductus symphaticus with the ductus cysticus. It is a reusable instrument intended for transient use.

Prostatic Dilators: A round, rod-like instrument for examining and dilating the urethra. It has a semi-rigid, conical shaft that merges into a narrow, bulbous finish. The instrument can be made hollow to be able to pass other, narrow instruments. It is a reusable instrument.

Vascular Dilators: A surgical instrument that is inserted into a vessel and used to dilate or adjust a vessel during a surgical procedure. The instrument is not intended for use on the blood vessels arteriae pulmonales, aorta ascendens, arcus aortae, descending aorta to bifurcation aortae, coronary arteries, common carotid artery, external carotid artery, internal carotid artery, cerebral arteries, brachiocephalic trunk, venae cordis, vena pulmonary, superior vena cava, and inferior vena cava. It is a rod-like instrument with a long, slender solid body and a distal tip in various shapes and sizes. It is available in a variety of sizes and degrees of softness and is typically made of metal. It is a reusable instrument intended for transient use.

Lachrymal Dilators: An ophthalmological instrument for use during ophthalmological procedures to dilate a cavity, canal or opening. Eye dilators also include instruments used for lacrimal or lacrimal canal surgery. It is a reusable instrument.

Needles:

Intestinal Needles; Muscular- & Skin suture; Reverdin Needles; Ligature Needles; post mortem needle: A cylindrically shaped, solid instrument of various diameters, whose distal tip can be made in various geometric shapes. It is used for introducing and pulling suture through superficial tissue, e.g. Skin or muscle, used. The suture serves to piece two wound edges together. It is a reusable instrument intended for transient use.

Eye Needles: A cylindrical, solid instrument of various diameters, whose distal tip can be made in various geometric shapes. It is used to insert and pull suture through ocular tissue. The seam connects two tissue surfaces. It is a reusable instrument intended for transient use.

Guide needles: A product for connecting a catheter to another product (e.g., a plug) or for introducing a flushing or drainage catheter into the body. It is a reusable instrument intended for transient use.

Spatula:

Tissue spatula: A surgical instrument, usually made of stainless steel, that is used to manipulate tissue or to remove material from a surface. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument intended for transient use.

Plaster Spatula: A instrument, usually made of stainless steel, for applying plaster. It is a reusable instrument

Abdominal-& Intestinal Spatulas; Organ spatulas: A surgical hand instrument for manipulating the intestine and abdominal organs during an abdominal surgical procedure. It typically has a large, flat, blunt blade (without sharp edges) at the distal end and a strong handle at the proximal end. It is usually made of stainless steel and is available in various designs, e.g. curved and flat, or straight and

flat with a rounded end, or right angled. Some models have a ridged leaf surface to better hold the slippery abdominal tissue. The handle is typically designed so that it can be kept non-slip. It is a reusable instrument intended for transient use.

Lung spatulas: A surgical instrument, usually made of stainless steel, which is used to manipulate lung tissue surfaces. The instrument consists of a proximal handle and a flat blade with no sharp edges at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. It is a reusable instrument intended for transient use.

Spatula: A surgical hand instrument, usually in spoon shape or rounded at the end of work. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. It is a reusable instrument intended for transient use.

Bladder Retractors; Gall Stone Scoops; Uterine Curettes: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon. It is a reusable instrument intended for transient use.

Tongue spatulas: A surgical hand-held tongue removal instrument to facilitate examination of the surrounding organs / tissues. It is typically a double-ended, flat product of various sizes with a slightly angled working end to better manipulate the tongue. It is made of stainless steel and is a reusable instrument.

Cement spatula: A hand instrument for applying orthopedic bone cement to a surface. It is usually made of stainless steel and consists of a proximal handle and a flat blunt blade (no sharp edges) at the distal end. The shaft or the blade of the instrument can be straight, curved or variously angled. The instrument is reusable.

Specula:

Ear Specula: A surgical instrument, usually composed of two articulated branches, available in a variety of sizes, shapes and contours. It has a blunt distal end and is used in insertion and opening to stretch or stretch tissue in the opening of the ear. It is used for examination or during procedures. It is usually made of stainless steel with a self-opening mechanism or an adjustable mechanism. It is a reusable instrument.

Mirrors:

Laryngeal Mirrors; Laryngeal mirror handles; Middle Ear Mirror: A surgical hand instrument with a polished surface, for reflecting sufficient undiffused light to obtain a visible image of an object during an ENT procedure. This mirror is mounted on a long, narrow handle, which is held by the doctor in the vicinity of the object to be observed. It is a reusable instrument intended for transient use.

Probes:

Rectal director: A surgical instrument used for exploring the tissues and structures of the rectum during a gastroenterological / urological surgical procedure. It is a reusable instrument intended for transient use.

Probes: A narrow, rod-like, flexible metal surgical instrument, usually with a blunt, beaded, or pointed tip, for examining fistulas, lumens, or wounds during an ENT procedure. The distal end, in front of the bead / tip, may be bent at right angles to the shaft. It is a reusable instrument intended for transient use.

Ligature Conductors: A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. Guides may also be instruments with ring handles, with straight or slightly curved working ends. The version with ring handles has a termination near the distal end. The working ends form short, grasping arms. It is a reusable instrument intended for transient use.

Gall Stone Probes: A surgical instrument consisting of a slender rod of flexible metal with a blunt buttoned tip used to examine (bite) the bile duct, e.g. is used on the union of the ductus sympathicus with the ductus cysticus. It is used during surgical procedures to examine or dilate bile duct narrowing. It is a reusable instrument intended for transient use.

Lachrymal Probes: A slender rod-shaped instrument made of flexible metal, with a blunt tip of tuber, used to examine eye tissue. It is a reusable product.

Others:

Mallets: A manual instrument for applying force (e.g., by beating) to another product (e.g., a chisel, pry bar, driver) during a procedure. It is available in different designs and materials. The shaft is usually made of metal (chromed iron or stainless steel), but may also be made of synthetic materials (e.g., tufnol). The hammer head is usually solid, annular, mostly double-ended, and is made of solid (iron / steel) or a shock-absorbing material (e.g., rubber, plastic, or tufnol). It is a reusable instrument.

Myoma Screw: A manual surgical hand instrument for detecting and holding a fibroid (i.e., a tumor of muscle tissue, usually in the uterus) during its surgical removal. It is typically a one-piece instrument, similar to a corkscrew with a pointed, spindled working end (for screwing into the myoma) and a long, slender shaft and a T-shaped handle. It is usually made of stainless steel. It is a reusable instrument intended for transient use.

Cotton Applicators: A narrow, rod-like, flexible metal instrument, usually with a blunt, beaded, or pointed tip, for examining fistulas, lumens, or wounds during an ENT procedure. The distal end, in front of the bead / tip, may be bent at right angles to the shaft. It is a reusable instrument intended for transient use.

Dilatationsbougies: A round, rod-like instrument for examining and dilating the urethra. It has a semi-rigid, conical shaft that merges into a narrow, bulbous finish. The instrument can be made hollow to pass other, narrow instrument can. It is a reusable instrument.

Table top cutters: A heavy duty hand instrument with specially designed powerful handles and a cutting jaw for cutting orthopedic wire, small diameter pins or cerclages. It has a scissor-type design with straight or curved handles and is made of stainless steel. It is available in various sizes and the cutting edges are usually short, strong and dull. The final part consists either of a single or double ratio screw connection to apply a greater force. The instrument, commonly referred to as a wire cutter, may be designed as a front or side cutter, some models have carbide inserts to extend shelf life. It is a reusable product.

Areola markers: A hand-operated surgical instrument used for removing blemishes or for sampling skin tissue. It usually consists of a handle with tips at the end of work in various configurations; some peaks include a closed cutting zone, edged by a cutting blade. This is a reusable instrument.

Breast Enlargement Sets: A surgical instrument with a sharp pyramidal or conical tip used to pierce body cavity. It is usually assembled and used with a compatible sleeve, filling its lumen, thereby allowing insertion of the instrument. After piercing, the trocar is withdrawn and provides a working channel into the body cavity. It is a reusable instrument.

Screwdriver: A tool with a shaft that fits into a screw head to tighten or remove this screw. The working end may e.g. for a slotted screw, Phillips screw, Phillips, Hexagon. The shaft end may have a handle for manual operation or a connection for a chuck for electric drive. It can have an indication of the torque. It is a reusable instrument intended for transient use.

Lead Hands: A personal protective equipment used to shield the entire hand of an operator or other from unnecessary primary and scattered radiation in the context of diagnostic and therapeutic procedures. Each finger is shielded individually. When used for X-ray and nuclear medicine, they usually consist of a liquid resistant outer covering surrounding a thin layer of lead or lead equivalent material. It is a reusable instrument.

Extension Bows: The product is part of a traction unit. It consists of metal bars with a curved frame, for attaching hooks or pins for exerting a pull on the patient. The traction unit exerts a tensile force to allow the body parts, e.g. the cervical or lumbar spine, by applying a strap around the head or pelvis.

The adjustable spreader bar adjusts the angle at which the belt exerts the tension. It is a reusable instrument.

Endarterektomiestripper; Varicose vein cutlery; Venous Instruments, Venous Hooks, Venous Extractors: A surgical instrument to cut a piece out of a vessel. A stripper can be used to excise part of a vein or artery except for the arterial arteries pulmonary, ascending aorta, aortic arc, descending aorta, aortic bifurcation, coronary artery, common carotid artery, external carotid artery, internal carotid artery, arteriae cerebrales, brachiocephalic trunk, venae cordis, vena pulmonales, superior vena cava, and inferior vena cava. Vessel Stripper is available in two different versions: 1) a flexible stainless steel cable with a peeling dome or disc at one end and a guide tip at the other end; and 2) a rigid rod terminating in a closed loop or loop, e.g., an external stripper. The instrument is reusable after proper cleaning has been done and is intended for transient use.

Measuring probes: A product / instrument used for comparative measurement in clinical use, e.g. inner and outer diameters, lengths, depths or thicknesses to measure. It is a non-calibrated product that is reusable and intended for transient use.

Ligature Guides; Wire strut guide: A surgical instrument for guiding wire or ligatures through tissue. Typical guides may have flexible or rigid working parts with hooks, clips or claws for holding the material while passing through the tissue. The working end can have a button or an eye. It is a reusable instrument intended for transient use.

Filling instrument: A instrument for applying filling and restorative material to the patient's mouth. It is a reusable product.

Sharp spoon: A spoon-shaped surgical instrument for use in gastroenterological and urological surgical procedures. This instrument group includes e.g. Biopsy spoon, ovule spoon, placenta spoon, uterine spoon, secretion spoon. It is a reusable instrument intended for transient use.

Ear Loops: A surgical ENT instrument, usually consisting of a canal or cannula and two fixed and sliding finger rings. A sling positions a circle of round, flexible wire around tissue in the ENT area and holds it in place. By moving the sliding finger ring you can control how much wire slides through the channel. The instrument is commonly used to exfoliate depleted foreign bodies. It is a reusable instrument intended for transient use.

3.2 Contraindications

No contraindications are known.

4 Warning

!	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
!	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
!	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
!	Please observe the additional information enclosed with the products.
!	Remove all protective sleeves and films prior to first using or preparation for use.
!	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
!	Avoid improper throwing or dropping of instruments
!	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
!	A visual inspection of the instrument for damage and contamination must be carried out before each use!
!	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
!	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
!	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
!	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do **not** use hard brushes (e.g. metal brushes or metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

The drying time during validation was 1 hour.

6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at **35-40 kHz, 5 minutes** at least.

V To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft materials are not usable for the ultrasound bath.

6.5 Manual cleaning

- !** Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes.
- Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

6.6 Mechanical cleaning

Based on international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.


Only completely demineralized water should be used for cleaning, neutralisation and rinsing, in accordance with the „Guidance Compiled by the DGKH (Germany Society for Hospital Hygiene), DGSV (German Society for Sterile Supply) and AKI (Working Group on Instrument Reprocessing) for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors“ (which refers to DIN EN ISO 15883-1 Point 6.4.2)

Flexible (complex) instruments with invisible surfaces must be pre-cleaned manually before mechanical cleaning.



We recommend for all push shafts, -and pipe shafts instruments and instruments whose surfaces are on top of each others during the cleaning (f.e. bone forceps and gouge forceps) a manually pre-cleaning for an optimal cleaning result without residues.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical. The chemical cleaning should take place at **40°C -60°C** for at **least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between **80 and 95°C**, with an **exposure time as outlined in EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

V Parameters used for the validation of preparation	
Pre-rinsing	1 minute with cold tap water
Cleaning	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
Neutralization	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
Post-rinsing	2 minutes with cold DI water
Disinfection	Temperature: 90 °C (A ₀ 3000)
	Soaking Time: 5 minutes

6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

V Drying was omitted in the validation (worst case condition)

7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement. Carry out a function check in accordance with the intended application of the instrument.

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

 Parameters used for the validation of steam sterilization	
Prevacuum	3 times
Sterilization temperature	132 °C
Sterilization time	1,5 minutes (half cycle method)
Drying time	20 minutes

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packagings for sterilized packaging.

 During validation the instruments were packaged in hospital common sterilization packagings (paper/film packagings) and steam sterilized.

9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

12 Manufacturer and servicing address

HEBUmedical GmbH
Badstraße 8
78532 Tuttlingen / Germany
Tel. +49 7461 94 71 - 0
Fax +49 7461 94 71 - 22
eMail: service@HEBUmedical.de
Web: www.HEBUmedical.de

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**Our vast range of products covers over 10 000 various instruments.
They can be found in Internet or request our catalog.**

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