# Einfache Schneidinstrumente / Simple cutting instruments

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**Operating Manual**  
Page 17

DEUTSCH

ENGLISH

## Inhalt

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
# 1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Validated Parameters
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Non sterile
	Keep away from sunlight
	Dry storage required
 Hinweis auf eIFU	(Electronic) instruction for use

# 2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them. By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

### 3 Scope

The use of simple cutting instruments is mainly in the following areas:

#### Knives

Plaster Knives : The instrument is used in various procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Amputating Knives: The instrument is used in operations for the amputation of a limb. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Micro-Instruments: The instrument is used for surgical procedures in gynecology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Eye Knives; Keratomes: The instrument is used for surgical procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Flap knife; Phalangeal knife; Interosseous knife; Post mortem knives; Resection knives; Autopsy knife; Micro Vessel Knives; Micro-laryngeal Instruments; Knives; Mucosa Knife; Knives, others: The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Cartilage knife; Periosteal knife: The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Rhinoplastic and Nasal Knives; Tonsil Knives, Dissectors, Enucleators, Retractors; Tympanum Knives; Trigeminal Knives; Ear knives, dissectors and raspatory; Adenoid currettes; Swivel Knives: The instrument is used during ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Dermal Punches: The instrument is used in plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

#### Scalpels

Operating Knives; Scalpel Handles: The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel

#### Chisels

Bone Chisels and Gouges; Nasal Osteotomes; flat chisel; Rhachiotom; Sternum chisel: The instrument is used in various surgical procedures in orthopedics. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Septum Gouges: The instrument is used during ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

#### Others

Adenotomes; Enucleators: The instrument is used during ENT procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Dermatom; Tendon Stripper: The instrument is used in various surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Circumcision Instruments: The instrument is used for urological surgical procedures. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Menisectomy Knives: The instrument is used during a menisectomy. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Foreign Body Instruments: The instrument is used for procedures in ophthalmology. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Separators: The instrument is used for various procedures in general or plastic surgery. The procedure may only be carried out by suitably trained and qualified specialist personnel.

Hand and foot instruments: The instrument is used for nail and foot care.

### 3.1 Intended Use

#### Knives

Amputating Knives: A heavy-duty surgical instrument with a handle and a single-sided or double-sided cutting blade of various shapes and sizes for use in limb amputations. It is a reusable instrument intended for transient use.

Flap knife; Phalangeal knife; Interosseous knife; Post mortem knives; Resection knives; Autopsy knife; Micro Vessel Knives; Micro-laryngeal Instruments; Knives; Mucosa Knife; Knives, others: A surgical instrument with a handle and a scalpel blade (not interchangeable) used by doctors to cut or prepare tissues with downward pressure. The instrument is made of stainless steel. It is a reusable instrument intended for transient use. For a scalpel with replaceable blade, see the corresponding generic product group.

Plaster Knives: A hand instrument with a short, flat, strong cutting blade with a sharp edge for cutting or shortening plaster casts. It is usually a one-piece product made of stainless steel. It is a reusable product.

Cartilage knife: A surgical instrument for cutting, scraping or shaping cartilage tissue. It is a reusable instrument.

Periosteal knife: A surgical instrument for cutting the periosteum. It is a reusable instrument.

Rhinoplastic and Nasal Knives; Tonsil Knives, Dissectors, Enucleators, Retractors; Tympanum Knives; Trigeminal Knives; Ear knives, dissectors and raspatory: A surgical hand instrument for cutting anatomical conditions during a general ENT procedure. It is made of stainless steel and has a handle in various designs, which merges into a narrow working end with a narrow, sharp blade. It is a reusable instrument intended for transient use.

Adenoid Curettes; Swivel Knives: A surgical hand instrument for cutting and excision of lymphoid tissue during an adenoidectomy. It is usually performed as a long, narrow instrument with a handle and a concave, spoon-like working part. The working part has a sharp edge for removing the tissue without injuring the surrounding muscles. It is usually made of stainless steel. It is a reusable instrument intended for transient use.

Eye Knives; Keratomes: An ophthalmic surgical instrument for placing precise cuts in ocular tissue during ophthalmic surgery on the eye and surrounding structures. It is usually designed as a one-piece instrument with a sharp, single-edged blade at the distal end and a handle at the proximal end. It can have different blade shapes, e.g. phaco (dull or sharp), crescent-shaped, straight, and micro-phaco (for cataract operations). It is usually made of stainless steel and some models may be equipped with diamond blades. It is a reusable instrument intended for transient use.

Micro Knives: A surgical instrument for cutting tissue that consists of a long slender metal handle that tapers towards the distal end. The internal radius of the hook has a cutting edge. It is a reusable instrument intended for transient use.

Micro-Instruments: A shaft-like surgical instrument whose distal end has the form of a hook or more parallel hooks for fixing or pulling uterine fibroids during a gynaecological procedure. It is typically a one-piece stainless-steel instrument with a handle in various designs at the proximal end. It is a reusable instrument intended for transient use.

**Dermal Punches:** A hand-operated surgical instrument used for removing blemishes or for sampling skin tissue. It usually consists of a handle with tips at the end of work in various configurations; some peaks include a closed cutting zone, edged by a cutting blade. This is a reusable instrument.

## Scalpels

**Operating Knives:** A surgical instrument with a handle and a scalpel blade (not interchangeable) used by doctors to cut or prepare tissues with downward pressure. The instrument is made of stainless steel. It is a reusable instrument intended for transient use.

**Scalpel Handles:** An interchangeable component of a scalpel that works as a handle and picks up a matching blade. Together, it provides a surgical instrument for cutting or dissecting tissues. The instrument can be made of stainless steel or brass. For one-piece scalpels, there is a suitable generic product group. It is a reusable instrument intended for transient use.

## Chisels

**Bone Chisels and Gouges; Nasal Osteotomes; flat chisel; Rhachiotom:** A surgical, chisel-like hand instrument for cutting and / or shaping bones by bumping during orthopedic surgery. It is held by the surgeon who applies a manual force to the proximal end of the instrument via a surgical hammer. The distal end (the cutting or sharp edge) is sharp, often flat, sometimes curved (concave), and usually folded on both sides. It is made of stainless steel. It is a reusable instrument intended for transient use.

**Sternum chisel:** A single-sided tapered blade surgical instrument for cutting and shaping bone during orthopedic surgery. The surgeon holds it in one hand while hitting the proximal end with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, flat or curved (concave), also called a hollow blade). It is made of high quality stainless steel and may consist of one piece, or have a synthetic handle (such as Tufnol). It is a reusable instrument intended for transient use.

**Septum Gouges:** A surgical instrument with a single sided beveled blade for cutting and shaping nasal bones during ear, nose and throat (ENT) operations. The surgeon holds it in one hand while beating the proximal end of the instrument with a surgical gavel in the other hand. The distal end (the cutting edge) is sharp, but often flat but may be curved (concave). It is typically one piece and made of high quality stainless steel. It is a reusable instrument intended for transient use.

## Others

**Dermatom:** A surgical hand instrument for cutting thin skin slices for transplantation or excision of small skin lesions. This requires a special cutting blade, which is used for this purpose in the instrument. It is a reusable instrument intended for transient use.

**Enucleators:** A surgical hand instrument for the dissection of the tonsils during an ENT procedure. It is usually spoon-like or rounded at the working ends. It is manufactured in various designs, e.g. with a shaft, which merges at the distal end in a rounded, bent or angled tip, or even double-ended. The working end may be tapered, flat, sharp or dull, straight or angled, serrated or smooth. It is made of stainless steel and is available in different sizes and shapes. It is a reusable instrument intended for transient use.

**Circumcision Instruments:** A surgical instrument for the controlled removal of the penis foreskin during circumcision. It is typically made of high quality stainless steel and may have various shapes, e.g. the bell clamp (bell-shaped clamp). The foreskin is pushed over the bell-shaped protective body and the glans are positioned therein. A screw mechanism is tightened so that the foreskin is compressed in an annular manner and fitted with e.g. a scalpel along the compression ring can be cut off. It is a reusable instrument intended for transient use.

**Adenotomes:** A surgical instrument designed to excrete hypertrophic lymphoid tissue in the nasopharynx (i.e., pharyngeal tonsils or polyps) during an adenoidectomy. It is usually a long, slender, manual hand instrument with cutting edges at the working end and a handle with a mechanism for

actuating the blades. It is available in various designs and sizes and can be a one-piece instrument or require the insertion of blades into the handle. It is a reusable instrument intended for transient use.

Meniscotomy Knives: A cutting, surgical instrument whose handle is T-shaped to be able to apply more force. The handle merges into a chisel-shaped edge. The instrument is used to cut crescent-shaped parts of the capsule or cartilage area. It is a reusable instrument intended for transient use.

Tendon Stripper: A surgical instrument consisting of a handle and a semi-circular, domed working end that ends in a sharp cutting edge. Tendon strippers are used to cut out a piece of ligament, tendon or fascia to transplant them. It is a reusable instrument intended for transient use.

Foreign Body Instruments: An ophthalmic instrument to remove foreign bodies from inside the eye or from the surface of the eyeball. It is typically made of high quality stainless steel and has a fine handle at the proximal end and a thin concave blade with a rounded nose at the distal end. The ophthalmologist uses it to remove the foreign body / object. It is a reusable instrument intended for transient use.

Separators: A surgical hand instrument, usually in the form of spoon or rounded at the working end. It is used to separate soft tissue or body structures in interventions of general or plastic surgery. It usually has a handle that continues into a shaft that has a tip at the distal end. The tip may be pointed or flat, sharp or dull, angled or straight on the shaft. The instrument is made of high quality stainless steel and available in various shapes and dimensions. The instrument is reusable.

Hand and foot instruments: A non-invasive instrument for podiatrists to treat deformities and diseases of the foot. It is a reusable instrument.

### 3.2 Contraindications

No contraindications are known.

## 4 Warning

	Medical products are delivered in a non-sterile condition and must be cleaned, disinfected and sterilized prior to their initial use.
	The use of faulty instruments is in principle forbidden and they have to go through the whole cleaning process before return.
	Please take into consideration that through higher power a bigger damage of the tissue can result: f.e. on forceps: the power at the end of the jaw is higher than at the tip of the jaw
	Please observe the additional information enclosed with the products.
	Remove all protective sleeves and films prior to first using or preparation for use.
	The safe combination of different products or of products with implants must be reviewed prior to clinical application by the user.
	Avoid improper throwing or dropping of instruments
	Avoid mechanical overstressing of the instrument beyond the structural design, this can lead to breakage and deformation!
	A visual inspection of the instrument for damage and contamination must be carried out before each use!
	To prevent all contact corrosion, instruments with damaged surfaces must be separated immediately.
	If the products are used on patients with transmissible spongiform encephalopathy or HIV infection, we decline any responsibility for their reuse.
	After ophthalmical use, please pay attention to water quality during treatment (according to the specifications of AAMI TIR34 and the recommendations of the Rober-Koch-Institute on preparation of medical devices)!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## 5 Handling

The type of treatment must be determined in each individual case by the surgeon in cooperation with the internist and the anaesthetist.

For operational use in various surgical disciplines must be done by appropriately trained and qualified personnel.

## 6 Preparation

The person in charge of preparatory treatment is responsible for ensuring that the treatment is duly carried out using the relevant equipment, materials and personnel in the treatment facility and so achieves the desired result. This necessitates validation and routine monitoring of the process used. We urge you to take note of the national regulations dealing with instrument preparation.

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

## 6.1 Reutilization restrictions

Frequent repeat preparatory treatment has minimal effects on the product. The end of the product life is normally determined by wear and damage due to use

## 6.2 Information on instrument preparation

- Use cleaning and/or disinfection agents with a pH-value within 9-10. Please observe manufacturer instructions regarding dosage, exposure time and renewal of solutions.
- Do **not** use hard brushes (e.g. metal brushes or metal sponges) or coarse abrasive cleaners.
- Never leave instruments in cleaning or disinfection agents for longer than the specified time.
- Only used demineralized water for rinsing.
- Rinse and dry carefully through channels and pipes.
- Sensitive instruments must be cleaned in a storage or clamping fixture.
- Observe manufacturer instructions of cleaning – and sterilizing equipment.

## 6.3 Preparation at the place of use

Directly after using remove coarse dirt of the instruments and rinse out the working cannulas. Do not use fixing agents or hot water (> 40°C), as this results in residues becoming fixed and can affect the success of the subsequent cleaning operation

Dismantle and/or open instruments as far as possible. Within short time after use the instruments clean the instruments for reducing a drying of the residues.

This enables an easier cleaning. If instruments come into contact with corroding medicines or cleaning agents, wash these up with water immediately after use.

Longer drying times, e.g. for dry disposal are not validated and therefore not recommended.

The drying time during validation was 1 hour.

## 6.4 Ultrasound bath (optional)

All instruments must be opened, dismantled and any cavities rinsed through.

Place instruments in the screen basket in such a way that overlaps and contact between instruments are avoided. Add cleaning agent to the water and adjust the temperature of the solution in line with the cleaning agent manufacturer's instructions.

The cleaning in the ultrasound bath should be at 35-40 kHz, 5 minutes at least.

To validate cleaning in an ultrasonic bath, the test items were ultrasonically treated in Neodisher Mediclean forte 0,5 % for 5 minutes.

Subsequently rinse instruments including all cavities before cleaning and disinfection.

Medicine products which are possessing a bad ultrasound transmission, e.g. soft materials are not usable for the ultrasound bath.

## 6.5 Manual cleaning

Since mechanical processes can be standardized, reproduced and therefore validated, mechanical cleaning/disinfection should be preferred to manual processes.

Manual cleaning and disinfection process is not validated and therefore needs to be validated additionally by the end user.

## 6.6 Mechanical cleaning

Based on international standards (EN ISO 15883) and national directives, only validated machine cleaning and disinfection methods may be used. For the mechanical cleaning we recommend a standard programme for surgical instruments, f.e. instruments from Miele.


Only completely demineralized water should be used for cleaning, neutralisation and rinsing, in accordance with the „Guidance Compiled by the DGKH (Germany Society for Hospital Hygiene), DGSV (German Society for Sterile Supply) and AKI (Working Group on Instrument Reprocessing) for the Validation and Routine Monitoring of Automated Cleaning and Thermal Disinfection Processes for Medical Devices as well as Advice on Selecting Washer-Disinfectors“ (which refers to DIN EN ISO 15883-1 Point 6.4.2)

Flexible (complex) instruments with invisible surfaces must be pre-cleaned manually before mechanical cleaning.



We recommend for all push shafts, -and pipe shaft instruments and instruments whose surfaces are on top of each others during the cleaning (f.e. bone forceps and gouge forceps) a manually pre-cleaning for an optimal cleaning result without residues.

Observe the following by loading:

- Place the dismantled/opened instruments securely in the tray.
- Instruments with openings and gaps have to be faced down with the opened side so that they can be cleaned and no water of the cleaning process is collecting inside them. If available use balanced devices for rinsing.

- Do not overload trays, avoid creating any overlaps.

Preliminary rinsing (cold, if applicable fully demineralized water without additives) is followed by chemical. The chemical cleaning should take place at **40°C -60°C** for at **least 5 minutes**.

We recommend products with a **pH-value within 9-10**, e.g. Neodisher MediClean forte from Dr. Weigert. The cleaning agents used should be selected depending on the material and properties of the instruments and in accordance to national regulations: If there is a high chloride concentration in the water, pitting and tension crack corrosion can occur on the instruments. The occurrence of this type of corrosion is minimized by using alkaline cleaning agents and demineralized water. By adding an acid-based neutralization medium, the rinsing off of alkaline cleaning agent residues is facilitated during the first intermediate rinsing process (warm or cold water).

In order to prevent the formation of deposits, it is advisable to use neutral cleaners where the water quality is unfavourable. After the second intermediate rinsing process, thermal disinfection takes place.

The thermal disinfection should take place at temperatures of between **80 and 95°C**, with an **exposure time as outlined in EN ISO 15883**.

After the finishing of the programme take the good out of the machine because corrosion can arise if the instrument remains in the machine.

 Parameters used for the validation of preparation	
<b>Pre-rinsing</b>	1 minute with cold tap water
<b>Cleaning</b>	Temperature: 55 °C
	Soaking Time: 5 minutes (worst case)
	Neodischer Mediclean forte 0,4% (worst case)
<b>Neutralization</b>	Temperature: cold DI water
	Soaking Time: 2 minutes
	Neodisher Z 0,1%
<b>Post-rinsing</b>	2 minutes with cold DI water
<b>Disinfection</b>	Temperature: 90 °C (A <sub>0</sub> 3000)
	Soaking Time: 5 minutes

## 6.7 Drying

Ensure adequate drying by the cleaning and disinfection device or using other suitable measures.

 Drying was omitted in the validation (worst case condition)

## 7 Maintenance, inspection

After cooling to room temperature, the instruments must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas. Instruments which are not free of residues must be returned for a complete retreatment process.

Carry out the functional check mentioned above. Instruments with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty instruments that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

### 7.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion.

If there are joints, the instruments should be oiled with a care product before the functional test. We recommend a medical white oil based on paraffin oil.

Check instruments with joints for ease of movement. Carry out a function check in accordance with the intended application of the instrument.

Essential tests for instruments with cutting edge are among others:

- Smooth cutting edges

Defective products must not be used and must have undergone the complete preparatory treatment process again before being returned.

## 8 Sterilization

Prior to sterilization, products must undergo cleaning and disinfection, be rinsed off without residue using demineralized water and subsequently dried. HEBUmedical recommends using a validated steam sterilization process (e.g. sterilizer in compliance with EN 285 and validated in accordance with DIN EN ISO 17665-1).

The validated parameters refer to reusable surgical instruments. The validated parameters should also be observed for the other products described, unless a different procedure is explicitly described.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

 <b>Parameters used for the validation of steam sterilization</b>	
<b>Prevacuum</b>	3 times
<b>Sterilization temperature</b>	132 °C
<b>Sterilization time</b>	1,5 minutes (half cycle method)
<b>Drying time</b>	20 minutes

The vapour must be free of ingredients, recommended limiting values of feed water and vapour condensate are determined through EN 285.

Other sterilization processes are compatible but not validated from HEBUmedical.

When loading, observe the recommended total weight. After the sterilization, check the sterile product packaging for damage, and inspect the sterilization indicators

## 8.1 Packaging

Compliant packaging of products for sterilization in line with ISO 11607. Packaging used must be suitable for the instruments and protect them from microbiological contamination during storage. The seal must not be under tension. HEBUmedical recommends container or hospital common sterilization paper/film packaging for sterilized packaging.

**V** During validation the instruments were packaged in hospital common sterilization packaging (paper/film packaging) and steam sterilized.

## 9 Lifetime

The steam sterilization procedure was validated by laboratory tests. The products were sterile validated at a pre-vacuum of at least 5min duration and a temperature of 134°C for a lifetime of 50 cycles.

You can continue to use the instruments at your own responsibility over this cycle value if the tests described in chapter 7 have been successfully completed.

## 10 Storage

Store products in a dry, clean and dust-free environment at moderate temperatures from 5°C to 40°C. Protect from the effects of the sun's rays and artificial light.

## 11 Warranty / Repair

Our products are manufactured from high-grade materials and carefully checked prior to dispatch. However, even if used properly in accordance with their intended purpose they are subject to a greater or lesser degree of wear depending on their intensity of use.

This wear is technically induced and unavoidable.

Should faults occur independently of wear, please contact our customer services. Defective products should no longer be used.

They must undergo the complete preparatory treatment process before being returned.

## 12 Manufacturer and servicing address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

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**Our vast range of products covers over 10 000 various instruments.  
They can be found in Internet or request our catalog.**

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HEBUmedical GmbH