HF 9945-xx
HF 9946-xx

HEBU HFiT[®]-T

Nadelelektroden / Needle Electrodes/ Elettrodi ad ago

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Seite 3

Operating Manual

Page 11

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Page 19

DEUTSCH

ENGLISH

ITALIANO


HF 9945-xx
HF 9946-xx

Inhalt

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HF 9945-xx
HF 9946-xx

CONTENT

1	Symbol descriptions	12
2	Introduction	13
3	Scope	13
3.1	<i>Intended Use</i>	13
3.2	<i>Contraindication</i>	13
3.3	<i>Maximum electrical capacity</i>	13
4	Instructions for Use	14
4.1	<i>Connection of the Needle Electrode</i>	14
5	Warnings	15
6	Safety Remarks	16
7	Preparation and Sterilization	17
8	Inspection	17
8.1	<i>Functional check</i>	17
9	Types of Needle Electrodes	18

1 Symbol descriptions

Symbol	Definition
	CE-labelling
	Attention
	Manufacturer
	Lot-description
	Reference code
	Medical device / FDA Prescription device
	Medical device
	Keep away from sunlight
	Dry storage required
	(Electronic) instruction for use
	No intervention on patients with cardiac pacemakers.
	Temperature
	sterile

Do not reuse

2 Introduction

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

By purchasing this instrument, you are now the owner of a high-quality product whose use and correct handling are described in the following. In order to minimize possible risks to patients and users, please observe these instructions carefully. Use, disinfection, cleaning and sterilization may only be performed by suitably trained specialist personnel.

3 Scope

Our products are exclusively intended for professional use by appropriately trained and qualified personnel and may only be acquired by them.

HEBU HFiT®-T Needle Electrode is mainly used in HF liver surgery for the treatment of inoperable tumors.

3.1 Intended Use

HFiT®-T needle electrodes serve as an applicator for radiofrequency thermal ablation of inoperable tumors of the liver and other pathogenic biological tissues. High-frequency current is introduced into the tissue via the HFiT®-T needle electrode and perfused with physiological saline solution. HFiT®-T needle electrodes must only be used with the HEBU HFiT®-T electrosurgery system. It is a disposable product.

3.2 Contraindication

The HF monopolar electrodes must not be used in conjunction with the central nervous or circulatory system. Patients with pacemakers and in vivo cardiac defibrillators may have restrictions on use. HF monopolar electrodes should not be used in the immediate vicinity of bone, tooth root and metal.

3.3 Maximum electrical capacity

Maximum electrical capacity: $0,5kV_p$.

4 Instructions for Use

- Before inserting HFiT[®]-T needle electrode into the target area, aeration of the infusion channel and flushing with physiological saline solution must be ensured, in order to prevent blood rising into the inside of the needle.
- Air bubble inclusions must be avoided.
- The insertion can take place under ultrasound, MRT or CT monitoring.
- If the need arises, any blockage caused by clotting inside the needle can be eliminated with the sterile obturator. If this is not successful, HFiT[®]-T needle electrode must be rejected.
- HFiT[®]-T needle electrodes are supplied sterile.

4.1 Connection of the Needle Electrode

- The HFiT[®]-T needle electrode is connected to a standard commercial syringe pump infusion line via the Luer-Lock connection.

The HFiT[®]-T needle electrode is contacted via the lateral socket with cable HF 9571-30

5 Warnings

The products may only be used in combination with HEBU accessories and only by clinically instructed and qualified personnel. Electrodes must not be twisted, as this may cause damage to the material and the electrode might break. Power output of the HF unit: Always select the lowest possible output of the HF unit for any procedure. If the coagulation capacity of the electrode is below normal, do not increase the output of the HF unit without a thorough preliminary check. In particular, the following must be checked: correct attachment of all HF cables and plugs, correct activation of current by the finger - or footswitches, undamaged insulation of the HF cables and instrument, and check that the distal end of the electrode is dean and undamaged. Before operating, check that the signal transmission of the monitors works without interference (e.g. without noise); this check is made by pressing the yellow button (CUT) and the blue button (COAG).

Always unplug the HF-cable by holding on the connector. **Never pull the cable itself**, as this may lead to defects. Cables with defective/brittle insulation or with a broken electrical wire (due to severe bending) might cause a fire. Active HF instruments and HF handles must not be placed on the patient in order to avoid injuries in case of involuntary HF activation or activated "Auto Start" mode.

Endogenous burns are those caused by high current density in the patient's tissue. Possible reasons are: The patient is inadvertently positioned in contact with electrically conductive parts. Direct contact between skin areas and HF cables may lead to electrical capacitance which in turn may cause burns.

Exogenous burns are those caused by heat of ignited fluids or gases. They may also be caused by explosions. Possible reasons include: Ignition of skin cleaning agents and disinfectants, ignition or narcotic gases etc.

Cardiac pacemakers may be damaged by electrosurgical current. Consult a cardiologist before operating. Never use electrosurgical current on outpatients fitted with pacemakers.

6 Safety Remarks

	HFiT [®] -T needle electrodes are intended for single use only!
	Do not use HFiT [®] -T needle electrodes when packaging is damaged!
	Do not load HFiT [®] -T needle electrodes by bending or twisting!
	Do not use device in explosive environments!
	Verify that combustible materials used in cleaning agents, disinfectants or solvents for adhesives have been completely removed prior to surgery.
	Please refer to the manual instruction of the HEBU HFiT [®] -T to obtain information about radio frequency thermo ablation.
	Observe manufacturer's instructions for electrical load capacity!
	Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

7 Preparation and Sterilization

HFiT®-T Needle electrodes are sterile packed and must not be cleaned and sterilized before use. They must be disposed of after single use.

Optionally, it is also possible to purchase the HFiT®-T needle electrodes non-sterile. The non-sterile HFiT®-T needle electrodes must be cleaned and sterilized before use as follows:

Rinse Needle electrode thoroughly with demineralized water.

On using the fractionated vacuum method, **sterilization** must be performed with at least **134°C (USA 132° C)** with a **minimum dwell period of 3 minutes**. Vacuum drying must then be carried out for at least 20 minutes.

Do not sterilize HFiT®-T Needle electrodes in hot air sterilizers.

This product is a disposable product and must **NOT** be reused!

8 Inspection

The products must be visually inspected for protein residues and other contamination, paying particular attention to cavities, blocks, inclusions, pipes, and other inaccessible areas.

Carry out the functional check mentioned above. Products with stains, which are blunt, bent, no longer function or which are otherwise damaged must be segregated!

To help identify faulty products that need to be sorted out, we recommend the brochure "Instrument Reprocessing" from the Working Group "Instrumenten Aufbereitung". This includes in particular Chapter 8 "Checks and Care" and Chapter 12 "Surface Changes: Deposits, Discoloration, Corrosion, Aging, Swelling and Stress Cracks".

8.1 Functional check

A newly purchased product must be subjected to a thorough visual and function check after its delivery and before each use. If a new product has externally visible defects or does not operate according to the instructions, contact the manufacturer or his responsible representative immediately.

When used as intended, these products are more or less subject to wear, depending on the intensity of use. This kind of wear is caused by technical factors and cannot be avoided. Therefore carry out a thorough visual examination and a function test prior to each application of these products. Carry out all plug connections only with the compatible articles and plug in the plug completely.

Products must be checked for irregularities. Paying attention to cracks, fractures and the occurrence of corrosion. Carry out a function check in accordance with the intended application of the product

Defective products must not be used!

9 Types of Needle Electrodes

Different types of HFiT®-T needle electrodes are available:

HEBU Item No.	Diameter	Length	Active Electrode
HF 9945-15	2,0 mm	150 mm	15 mm
HF 9945-20	2,0 mm	150 mm	20 mm
HF 9946-15	2,0 mm	200 mm	15 mm
HF 9946-20	2,0 mm	200 mm	20 mm
HF 9946-25	2,0 mm	200 mm	25 mm
HF 9946-30	2,0 mm	200 mm	mm

10 Storage

Storage of the sterilized instruments in a dry, clean and dust-free environment. Products must be protected from direct sunlight in storage.

11 Warranty/Repair

Our products are made of high-quality materials and are carefully checked before delivery. However, even when used as intended, they are subject to wear to a greater or lesser extent depending on the intensity of use.

This wear is due to technical reasons and is unavoidable.

However, if faults occur that are not related to wear, please contact our customer service department. Defective products may no longer be used and must have undergone the entire preparatory process before being returned.

12 Servicing and manufacturer address

Should you require the instructions for use in paper form, please use the contact details below. The instructions for use in paper form will be made available to you within seven calendar days of receipt of the request.

Alternatively, you can print out the electronic instructions for use yourself.

HEBUmedical GmbH
 Badstraße 8
 78532 Tuttlingen / Germany
 Tel. +49 7461 94 71 - 0
 Fax +49 7461 94 71 - 22
 eMail: service@HEBUmedical.de
 Web: www.HEBUmedical.de

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HEBUMedical GmbH